MedixMAB"

Technical Note

medixbiochemica.com



Contents

CRP	3
SAA	6
Calprotectin	7
PCT	9
NGAL	12

Introduction

Inflammation is part of the immune system's natural response to infection or tissue injury. Acute inflammation serves as a mechanism of initiating the elimination of infectious agents or damaged tissue by recruiting blood cells to kill invading micro-organisms, which in turn promotes wound healing. Classical symptoms of acute inflammation are redness of the skin, swelling, heat, and pain. If the inflammation persists and becomes chronic, it can eventually cause chronic diseases including atherosclerosis, cancer, Alzheimer's or autoimmune diseases. For this reason, inflammation has started to become one of the most studied areas of medical research.

One of the most common biomarkers for diagnosing inflammation is C-reactive protein (CRP), which is measured from blood samples and indicates whether the inflammation is caused by bacteria. If CRP values rise, doctors will prescribe antibiotics to treat a bacterial infection. If the inflammation is caused by viral infection, CRP levels do not normally increase and therefore antibiotics are not needed. High-sensitivity CRP (hs-CRP) is sometimes used to measure chronic inflammation as it is more precise than standard CRP when measuring baseline (i.e., normal) concentrations of CRP. Recently, other markers such as procalcitonin (PCT) and serum amyloid A (SAA) have also become more widely used for diagnosis of acute inflammation.

Medix Biochemica has concentrated on inflammation biomarkers as a key focus area. For over 35 years, we have continually been at the forefront of monoclonal antibody development, using modern *in vitro* production methods with excellent batch-to-batch consistency.

CRP

C-Reactive protein (CRP) has been named according to its capacity to precipitate the C-polysaccharide of *Streptococcus pneumoniae.* CRP is secreted in response to a variety of inflammatory cytokines, and its main biological role is the activation of the complement system and of other proinflammatory processes. More specifically, CRP binds to phosphocholine and related substances, which are widely distributed in polysaccharides of pathogens and in cellular membranes. CRP therefore recognizes a range of pathogenic targets as well as membranes of damaged and dying host cells, making it one of the key activators of the classical complement pathway.¹⁴

The levels of CRP in blood increase rapidly in response to trauma, inflammation, and infection. It is thus not a specific marker for a particular disease, but rather a widely used indicator of bacterial or fungal infection, or various inflammatory conditions, such as rheumatoid arthritis or inflammatory bowel disease. Normal levels of CRP in blood are generally less than 3 mg/L. In the case of acute inflammation or bacterial infection, CRP levels can rise dramatically, up to 10,000-fold. In recent years, studies have indicated that slightly elevated CRP levels might be linked to several other disease types as well, such as cardiovascular diseases or cancer. However, the physiological connection between CRP levels and the risk of cardiovascular events or malignant tumors is not yet completely understood.⁵⁻⁶

The human CRP molecule is a pentamer, i.e. it is composed of five identical polypeptide subunits, each containing 206


amino acid residues and having a molecular mass of 23 kDa. It belongs to pentraxins, which are an evolutionary conserved family of proteins with a unique architecture – a planar, doughnut-shaped ring.⁷

Medix Biochemica's extensive experience on monoclonal antibodies used in CRP diagnostics extends back to the 1980's when our first products – Anti-human CRP mouse monoclonal antibodies 6404 and 6405 – were developed. Since then, the product family has been extended, and it currently includes five different mouse monoclonal antibodies recognizing CRP.

The CRP antibodies have been studied using multiplexed surface plasmon resonance (SPR) technology and ELISA to identify good antibody pairs for CRP detection. It should be noted that the functionality of two antibodies as a pair is dependent on the method and conditions used in the assay.

The association and dissociation characteristics of an antibody-antigen binding reaction can have significant effects on the usability of antibodies in IVD applications. The association rate constant (on-rate) and the dissociation rate constant (off-rate) for Medix Biochemica's anti-CRP antibodies have been measured using SPR. Anti-human CRP antibodies have no cross-reactivity with the serum amyloid P (SAP) component.

Scientific publications: page 14



The human CRP molecule is a pentamer.

CRP antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
6402	100145	5	24	lgG ₁	ELISA, IT, LF
6403	100146	5	36	IgG ₁	ELISA, IT, LF
6404	100058	2.0-2.5	24	lgG ₁	ELISA, IT, LF, WB
6404	100061	> 6	24	lgG ₁	ELISA, IT, LF, WB
6405	100358	5	18	IgG ₁	ELISA, IT, LF, WB
6407	100147	5	24	lgG ₁	ELISA, IT, LF

Anti-human CRP monoclonal antibodies


Pair recommendations

				Detection		
		6402	6403	6404	6405	6407
	6402	+	+	-	+	+
ย	6403	+	+	-	+	-
Capture	6404	-	-	+	+	-
Ŭ	6405	+	-	-	+	-
	6407	+	-	-	+	+

Kinetic parameters

CRP antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _p (M)
6402	7.3 x 10⁵	1.2 x 10 ⁻⁴	1.6 x 10 ⁻¹⁰ = 0.16 nM
6403	9.6 x 10⁵	1.6 x 10 ⁻⁴	1.6 x 10 ⁻¹⁰ = 0.16 nM
6404	7.7 x 10 ⁶	4.1 × 10 ⁻⁴	5.3 x 10 ⁻¹¹ = 0.05 nM
6405	2.6 x 10 ⁶	6.6 x 10 ⁻⁵	2.6 x 10 ⁻¹¹ = 0.03 nM
6407	1.1 x 10 ⁶	1.3 x 10 ⁻⁴	1.1 x 10 ⁻¹⁰ = 0.11 nM

Inflammation



CRP standard curves

5



Acute-phase serum amyloid A (SAA) proteins are synthesized in response to inflammatory signals, similar to CRP. SAA1 and SAA2 – the two major types of SAA found in humans – are highly homologous proteins with a molecular mass of about 12 kDa. When released into circulation, SAAs associate with high-density lipoprotein (HDL) particles. Several immunological roles have been suggested for SAAs, such as the induction of cytokine synthesis and recruitment of immune cells to inflammatory sites, however the precise physiological function of SAAs is still not fully understood.⁸⁴⁰

SAA concentrations in blood increase within a few hours after the onset of inflammation. The normal levels vary from <2–8 mg/L, and the concentration can increase by as much as 1000-fold during acute inflammation. Like CRP, SAA can therefore be used as an indicator for microbial infections or for various inflammatory conditions. An increasing amount of evidence has also been collected about the significance of chronic inflammation in the development and progression of many other severe diseases, such as various cancers and cardiovascular diseases, and SAA has been suggested to be a good prognostic marker for these patients.¹¹⁻¹³ In addition, as SAA proteins are highly conserved in vertebrates, SAA measurements are also used in veterinary diagnostics. Medix Biochemica offers two mouse monoclonal antibodies for SAA detection, which can be used as a pair in diagnostic tests.

In addition to antibodies, Medix Biochemica offers a recombinant human SAA protein antigen, product code 610070, which contains the amino acids 19–122 (according to UniProt #P0DJI8; no signal peptide) and an N-terminal histidine tag.

Purity	97 kDa –	
Denaturing SDS- PAGE analysis of	66 kDa –	
recombinant human SAA protein antigen.	45 kDa –	
	30 kDa –	
	20.1 kDa –	
	14.4 kDa-	

Anti-human SAA monoclonal antibodies and recombinant antigen

SAA antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
2201	100279	5	18	lgG _{2a}	ELISA, LF, IT
2203	100289	5	18	lgG ₁	ELISA, LF, IT
SAA antigen		Product code			
Descent to a b	CA 4 400	(10070			

Recombinant SAA, 100 µg 610070

Kinetic parameters

SAA antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _D (M)
2201	1.9 x 10⁵	1.8 x 10 ⁻⁴	9.6 x 10 ⁻¹⁰ = 1.0 nM
2203	1.1 x 10⁵	3.0 x 10 ⁻⁴	2.8 x 10 ⁻⁹ = 2.8 nM

Calprotectin

Calprotectin is a calcium-binding heterodimer formed by two proteins, S100A8 and S100A9, associating together. It has a molecular mass of 24 kDa, and each subunit binds two calcium ions. In the presence of abundant calcium, two calprotectins may also bind together to form heterotetramers.¹⁹

Calprotectin is secreted by neutrophils at the site of inflammation. The protein has been shown to play a role in regulating cytokine synthesis, and also possesses antimicrobial activity mediated through the sequestration of essential nutrient metals such as zinc and manganese.²⁰⁻²¹

Calprotectin levels in plasma and stool are significantly increased in infectious and inflammatory conditions. Analysis

of fecal calprotectin is commonly used to diagnose intestinal inflammations, especially inflammatory bowel disease (IBD). Information about fecal calprotectin levels provides valuable support for treatment decisions as it can consistently differentiate IBD from irritable bowel syndrome.²²

Medix Biochemica offers five different monoclonal antibodies for sensitive calprotectin detection, with specificities for both S100A8/9 subunits. In addition, Medix Biochemica's product portfolio includes three types of recombinant calprotectin antigens.

Anti-human calprotectin monoclonal antibodies and recombinant antigens

Calprotectin antibody	Product code	Calprotectin subunit binding specificity	Shelf life (months at +2–8°C)	Subclass	Applications tested
3403	100460	S100A9	12	IgG _{2a}	ELISA, LF
3404	100468	S100A8	24	IgG ₁	ELISA, LF
3405	100469	S100A8	N/D	IgG ₁	ELISA, LF
3406	100470	S100A8	24	IgG ₁	ELISA, LF
3407	100618	S100A8/S100A9 complex	18	lgG ₁	ELISA, LF

Calprotectin antigens	Product code	Subunit structure
Recombinant Calprotectin, 100 µg	610061	Heterodimer S100A8/9
Recombinant Calprotectin S100A8, 50 µg*	710018	S100A8 subunit
Recombinant Calprotectin S100A9, 50 µg*	710019	S100A9 subunit

* Also available in 500 µg and 1000 µg package sizes.


Pair recommendations

				Detection	l	
		3403	3404	3405	3406	3407
	3403	-	+	+	+	-
e	3404	+	-	-	-	+
Capture	3405	+	-	-	-	+
ü	3406	+	-	-	-	+
	3407	-	+	+	+	+

Kinetic parameters

Calprotectin antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _p (M)
3403	1.9 x 10 ⁶	9.1 x 10 ⁻⁴	4.8 x 10 ⁻¹⁰ = 0.5 nM
3404	4.8 x 10 ⁵	6.1 x 10 ⁻⁵	1.3 x 10 ⁻¹⁰ = 0.1 nM
3405	7.4 x 10⁵	1.5 x 10 ⁻⁴	2.0 x 10 ⁻¹⁰ = 0.2 nM
3406	7.6 x 10⁵	1.1 × 10 ⁻⁴	1.4 x 10 ⁻¹⁰ = 0.1 nM
3407	4.0 x 10 ⁴	5.0 x 10 ⁻⁵	1.3 x 10 ⁻⁹ = 1.3 nM

Binding properties



• 3403 recognizes both the S100A9 subunit and the S100A8/S100A9 complex

• 3404, 3405 and 3406 recognize the S100A8 subunit and the S100A8/S100A9 complex

• 3407 only recognizes the S100A8/S100A9 complex but not the individual subunits.

PCT

Procalcitonin (PCT) is a precursor protein of the hormone calcitonin, produced by the thyroid gland. It is a 116 amino acid protein that is cleaved intracellularly to form three peptides: N-terminal procalcitonin, calcitonin, and katacalcin. In healthy individuals, only the calcitonin polypeptide is secreted into the bloodstream, and PCT serum levels are very low. However, in patients with systemic infections, and especially sepsis, PCT is secreted by several cell types, not only in the thyroid gland, and levels rise rapidly. High PCT levels have also been identified in patients suffering from cardiogenic shock, severe systemic inflammatory response syndrome (SIRS), and trauma.¹⁴⁻¹⁶

The immunological roles of PCT have been shown to be linked to the neutralization of bacterial lipopolysaccharides (LPS) and to the reduction of cytokine release. The LPS are strong endotoxins and mediators of septic shock, and it has been suggested that PCT may neutralize LPS via the polycationic motifs present in its molecular structure.¹⁷

The association of elevated serum PCT levels with bacterial infections makes PCT a very promising biomarker for sepsis, as it can be used to discriminate bacterial from noninfective causes of systemic inflammations. Additionally, it has been suggested that PCT measurement is useful in the management of antibiotic therapy, which is becoming increasingly important due to the spread of antibiotic-resistant pathogens.^{15,16,18}

Medix Biochemica has several monoclonal antibodies for PCT, against both katacalcin and calcitonin subunits. We also offer the full-length PCT antigen.

Purity	97 kDa –	
Denaturing SDS-PAGE analysis of recombinant human PCT protein antigen.	66 kDa –	
	45 kDa –	
	30 kDa –	
	20.1 kDa –	d'
	14.4 kDa –	-

Anti-human PCT monoclonal antibodies and recombinant antigen

PCT antibodies	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Specificity	Applications tested
4003	100562	5	N/D	lgG ₁	Katacalcin	ELISA, LF
4004	100563	5	24	IgG ₁	Calcitonin	ELISA, LF
4005	100564	5	24	IgG ₁	Calcitonin	ELISA, LF
4006	100567	5	N/D	IgG ₁	Katacalcin	ELISA, LF
4008	100769	5	N/D	IgG ₃	Katacalcin	ELISA, LF
4051	700020	>1	60	IgG ₁	Calcitonin	ELISA, LF

PCT antigen	Product code
Recombinant PCT, 100 µg	610080

Pair recommendations

	Detection					
	4003	4004	4005	4006	4008	4051
4003	-	+	+	_	-	+
4004	+	-	-	+	+	-
4005	+	-	-	+	+	-
4006	-	+	+	-	-	+
4008	-	+	+	-	-	+
4051	+	-	-	+	+	-
	4004 4005 4006 4008	4003 - 4004 + 4005 + 4006 - 4008 -	4003 - + 4004 + - 4005 + - 4006 - + 4008 - +	4003 4004 4005 4003 - + + 4004 + - - 4005 + - - 4006 - + + 4008 - + +	4003 4004 4005 4006 4003 - + + - 4004 + - - + 4005 + - - + 4006 - + - - 4006 - + - - 4006 - + + -	4003 4004 4005 4006 4008 4003 - + + - - 4004 + - - + + + 4005 + - - + + + 4006 - + - - - +

Kinetic parameters

	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _D (M)
4003	4.8 x 10 ⁴	2.2 x 10 ⁻⁴	4.6 x 10 ⁻⁹
4004*	2.5 x 10 ⁴	-	-
4005*	2.9 x 10⁴	-	-
4006	6.7 x 10 ⁴	4.0 x 10 ⁻⁴	5.2 x 10 ⁻¹⁰
4008	N/D	N/D	N/D
4051	N/D	N/D	N/D


*The antibody does not dissociate in the reaction conditions used.

PCT binding properties

N-terminal part	Calcitonin	Katacalcin	
	4004	4003	4006
	4005		4008
	4051		

Inflammation

PCT standard curve



11

NGAL

Neutrophil gelatinase-associated lipocalin (NGAL) is a small glycoprotein belonging to the superfamily of lipocalins. NGAL is expressed in epithelial tissues associated with anti-microbial defense, and in the distal tubules and collecting duct of the kidney. Elevated levels of blood and urine NGAL are associated with acute or chronic renal failure, and elevated serum levels with active inflammatory bowel disease (IBD) and cardiovascular events. Medix Biochemica offers four mouse monoclonal antibodies for NGAL detection as well as a recombinant NGAL antigen. The protein is produced in *Escherichia coli* host cells and contains the amino acids 21–198 (according to UniProt #P80188) and a C-terminal histidine tag.

Anti-human NGAL monoclonal antibodies and recombinant antigen

NGAL antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
4202	100579	5	24	lgG ₁	ELISA, IT
4203	100580	5	24	lgG ₁	ELISA, IT
4204	100581	5	36	IgG ₁	ELISA, IT
4205	100582	5	24	lgG ₁	ELISA, IT


NGAL antigen	Product code
Recombinant NGAL, 100 µg	610012

Pair recommendations

		Detection			
		4202	4203	4204	4205
	4202	-	+	+	+
Capture	4203	+	-	+	+
Capl	4204	+	+	-	+
	4205	+	+	+	-

Purity


Denaturing SDS-PAGE analysis of recombinant human NGAL protein antigen.



NGAL antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _p (M)
4202	6.4 x 10 ⁵	4.0 x 10 ⁻⁴	6.3 x 10 ⁻¹⁰
4203*	2.2 x 10 ⁵	-	-
4204	1.8 x 10 ⁵	3.4 x 10 ⁻⁵	1.8 x 10 ⁻¹⁰
4205	5.6 x 10⁴	4.9 x 10 ⁻⁵	8.6 x 10 ⁻¹⁰

Kinetic parameters

*The antibody does not dissociate in the reaction conditions used.



NGAL standard curves

Page 3

Scientific articles describing the use of MedixMAB CRP antibodies include:

- Käpyaho K, Welin MG, Tanner P, Kärkkäinen T and Weber T (1989) Rapid determination of C-reactive protein by enzyme immunoassay using two monoclonal antibodies. Scand J Clin Lab Invest. 49:389-393.
- Elschner T, Scholz F, Miethe P and Heinze T (2014) Rapid Flow Through Immunoassay for CRP Determination Based on Polyethylene Filters Modified with ω–Aminocellulose Carbamate. Macromol Biosci. 14:1539-1546.
- Tarkkinen P, Palenius T and Lövgren T (2002) Ultrarapid, ultrasensitive one-step kinetic immunoassay for C-reactive protein (CRP) in whole blood samples: Measurement of the entire CRP concentration range with a single sample dilution. Clin Chem. 48:269-277.

References:

- Du Clos TW (2000) Function of C-reactive protein. Ann Med. 32:274–278.
- Krüger S & Welte T (2012) Biomarkers in Community-acquired Pneumonia. Expert Rev Resp Med. 6:203–214.
- 3. Pepys MB & Hirschfield GM (2003) C-reactive protein: a critical update. J Clin Invest. 111:1805–1812.
- Venugopal SK, Devaraj S & Jialal I (2005) Macrophage Conditioned Medium Induces the Expression of C-Reactive Protein in Human Aortic Endothelial Cells : Potential for Paracrine/Autocrine Effects. Am J Pathol. 166:1265–1271.
- The Emerging Risk Factors Collaboration (2012) C-Reactive Protein, Fibrinogen, and Cardiovascular Disease Prediction. New Engl J Med. 367:1310–1320.
- 6. Allin KH, Nordestgaard BG, Zacho J, et al. (2010) C-Reactive Protein and the Risk of Cancer: A Mendelian Randomization Study. JNCI J Natl Cancer Inst. 102:202–206.
- Du Clos TW (2013) Pentraxins: Structure, Function, and Role in Inflammation. ISRN Inflammation 2013, Article ID 379040, 22 pages.
- Uhlar CM & Whitehead AS (1999) Serum amyloid A, the major vertebrate acute-phase reactant. Eur J Biochem. 265:501–523.
- Eklund KK, Niemi K, Kovanen PT (2012) Immune functions of serum amyloid A. Crit Rev Immunol. 32:335–348.
- Christenson K, Björkman L, Ahlin S, et al. (2013) Endogenous Acute Phase Serum Amyloid A Lacks Pro-Inflammatory Activity, Contrasting the Two Recombinant Variants That Activate Human Neutrophils through Different Receptors. Front Immunol. 4:92. doi:10.3389/fimmu.2013.00092.
- Sung HJ, Ahn JM, Yoon YH, et al. (2011) Identification and validation of SAA as a potential lung cancer biomarker and its involvement in metastatic pathogenesis of lung cancer. J Proteome Res. 10:1383–1395.
- Meng YQ, Cao X, Wen ZS, et al. (2014) Preoperative level of serum amyloid A is superior to C-reactive protein in the prognosis of esophageal squamous cell carcinoma. Dis Esophagus 27:670–677.

- Zewinger S, Drechsler C, Kleber ME, et al. (2015) Serum amyloid A: high-density lipoproteins interaction and cardiovascular risk. Eur Heart J. 36:3007–3016.
- Maruna P, Nedelníková K & Gürlich R (2000) Physiology and genetics of procalcitonin. Physiol Res. 49 Suppl 1:S57–61.
- Hentschel D, Olivier CB, Bode C, Diehl P (2014) The Role of Procalcitonin in Septic Patients – A Brief Overview. Transl Med. 4:144. doi: 10.4172/2161–1025.1000144.
- Wacker C, Prkno A, Brunkhorst FM & Schlattmann P (2013) Procalcitonin as a diagnostic marker for sepsis: a systematic review and meta-analysis. The Lancet Inf Dis. 13:426–435.
- Matera G, Quirino A, Giancotti A, et al. (2012) Procalcitonin neutralizes bacterial LPS and reduces LPS-induced cytokine release in human peripheral blood mononuclear cells. BMC Microbiol. 12:68. doi: 10.1186/1471–2180–12–68.
- Quenot JP, Luyt CE, Roche N, et al. (2013) Role of biomarkers in the management of antibiotic therapy: an expert panel review II: clinical use of biomarkers for initiation or discontinuation of antibiotic therapy. Ann Intensive Care. 3(1):21. doi: 10.1186/2110–5820–3–21.
- Stríz I & Trebichavský I (2004) Calprotectin a pleiotropic molecule in acute and chronic inflammation. Physiol Res. 53:245–253.
- Liu JZ, Jellbauer S, Poe AJ, et al. (2012) Zinc Sequestration by the Neutrophil Protein Calprotectin Enhances Salmonella Growth in the Inflamed Gut. Cell Host & Microbe 11:227–239.
- Dorosz SA, Ginolhac A, Kähne T, et al. (2015) Role of Calprotectin as a Modulator of the IL27-Mediated Proinflammatory Effect on Endothelial Cells. Mediators of Inflammation, Article ID 737310, 16 pages.
- 22. Konikoff MR & Denson LA (2006) Role of fecal calprotectin as a biomarker of intestinal inflammation in inflammatory bowel disease. Inflamm Bowel Dis. 12:524–534.

More information on our products **medixbiochemica.com**



Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com

CPS = Counts per second CLIA = Chemiluminescence immunoassay ELISA = Enzyme-linked immunosorbent assay IT = Immunoturbidimetry LF = Lateral flow N/A = Not Applicable N/D = Not Determined WB = Western blot

Copyright © 6/2020 Medix Biochemica. All rights reserved.

MedixMAB and MedixAntigens are trademarks of Medix Biochemica. Medix Biochemica reserves the right to make changes and improvements to any of the products described in this document without prior notice.