

Technical Note

medixbiochemica.com



Contents

TSH	3
hCG	5
LH	8
FSH	10
Alpha subunit	11
SHBG	12
Progesterone	13
Prolactin	14
T4	15
Thyroglobulin	16
Growth hormone	17

Introduction

Hormones are a diverse group of signaling molecules produced mainly by the endocrine glands. They control a vast array of essential processes ranging from reproduction, sexual differentiation, development and growth to the regulation of metabolism and maintenance of cellular homeostasis. Due to the large diversity of cellular and bodily processes controlled by hormones, there are also various disorders related to defective endocrine functions. For example, several tumors, including thyroid, testicular and ovarian tumors, are assumed to be endocrine-dependent and have the potential to promote their own growth by excess hormone production.

Clinical manifestations of endocrine disorders, however, are often indefinite and nonspecific. Accordingly, laboratory measurements are usually required to diagnose disease and to monitor therapeutic outcomes. Since hormone levels in bodily fluids are relatively low, reliable quantification of hormones often requires sensitive assays, which are usually immunoassays based on monoclonal antibodies that specifically recognize the target analyte. In addition to disorders and malignancies, hormone levels can also be used to monitor normal bodily processes. For example, immunoassays for human chorionic gonadotropin (hCG) are used for pregnancy detection and monitoring, while different hCG forms can also be used to detect trophoblastic tumors or pregnancy-related disorders. Human luteinizing hormone (LH) is often used for predicting ovulation, but it can also be used to investigate menstrual irregularities or to diagnose puberty-related disorders.

Medix Biochemica has over 30 years of experience in producing premium-quality monoclonal antibodies for detection of hCG, LH, and other hormones. The company's optimized, industrial-scale *in vitro* production methods, certified batch-to-batch consistency, as well as expert customer service have made Medix Biochemica one of the most important antibody suppliers for the IVD community.

TSH

Thyroid-stimulating hormone (TSH, also known as thyrotropin) is a glycoprotein hormone secreted by the anterior pituitary gland. In response to TSH, the thyroid gland secretes the hormone thyroxine (T4) that is converted into triiodothyronine (T3) in the liver and other organs. As T3 affects several cellular metabolic pathways, TSH is considered a key regulator of the normal development and metabolism of the body.

TSH is an approximately 30 kDa heterodimeric protein formed by two noncovalently linked subunits, alpha and beta. The alpha subunit is shared with human FSH, hCG, and LH. The structure is stabilized by a unique segment of the beta subunit, called the seatbelt, which is wrapped around the alpha subunit.^{1,2}

In clinical diagnostics, measurements of serum TSH are most often used to diagnose and manage thyroid disorders, including both hyperthyroidism (overactive thyroid) and hypothyroidism (underactive thyroid). Reported TSH levels have varied historically depending on the test method or platform used, and different organizations have recommended different reference intervals for normal TSH levels. In general, levels between 0.4 to 4 mIU/L are now considered normal in adults, excluding pregnant women. Recently, the International Federation of Clinical Chemistry and Laboratory Medicine (IFCC) has initiated standardization and harmonization efforts, with very promising results.³

Medix Biochemica has been a leading provider of monoclonal antibodies recognizing human TSH for more than 30 years. We offer a wide selection of mouse monoclonal antibodies against both the alpha subunit and the full human TSH.


An independent laboratory conducted a large TSH study using 40 human clinical samples and all commercially available anti-TSH antibodies. Of all the tested antibodies, the best sensitivity was achieved with Medix Biochemica's TSH monoclonal antibodies (capture antibody 5405 or 5407 and label antibody: 5409, unpublished data).

Scientific publications: page 18

antibody co	ode		Shelf life (months at +2–8°C)	Subclass	Applications tested
5401 10	00023	1	12	lgG ₁	ELISA,CLIA
5404 10	00026	5	36	lgG ₁	ELISA, CLIA, LF
5405 10	00819	5	36	lgG _{2a}	ELISA, CLIA, LF
5407 10	00254	5	24	lgG ₁	ELISA, CLIA, LF
5408 10	00033	1	36	lgG ₁	ELISA, CLIA, LF
5409 10	00034	5	36	lgG ₁	ELISA, CLIA, LF

Pair recommendations

		Detection							
		5401	5404	5405	5407	5408	5409*		
	5401	-	-	+	+	-	+		
	5404	-	-	+	+	-	+		
ture	5405	-	-	-	-	+	+^A		
Capture	5407	-	-	-	-	+	+^A		
Ŭ	5408	-	-	+	+	-	+		
	5409*	-	-	+	+	_	_		




* Recognizes an epitope at the junction of alpha and beta subunits in hTSH A) Recommendation for detection of TSH with highest sensitivity

Kinetic parameters

TSH antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _p (M)
5401	3.6 x 10 ⁵	7.0 x 10 ⁻⁵	2.0 x 10 ⁻¹⁰ = 0.20 nM
5404	8.5 x 10 ⁵	3.9 x 10 ⁻⁵	4.6 x 10 ⁻¹¹ = 0.05 nM
5405	1.2 x 10 ⁶	1.9 x 10 ⁻⁴	1.5 x 10 ⁻¹⁰ = 0.15 nM
5407	1.1 x 10 ⁶	2.0 x 10 ⁻⁴	1.8 x 10 ⁻¹⁰ = 0.18 nM
5408	9.6 x 10 ⁵	3.8 x 10 ⁻⁵	4.0 x 10 ⁻¹¹ = 0.04 nM
5409	2.1 x 10 ⁶	1.0 × 10 ⁻⁴	4.8 x 10 ⁻¹¹ = 0.05 nM

Epitope regions

TSH antibodies can be divided into different epitope groups. Antibodies in the same group detect the same or overlapping epitopes.



hCG

Human chorionic gonadotropin (hCG) is a glycoprotein hormone essential for the maintenance of pregnancy. Physiologically, this placental hormone is produced and secreted mainly by embryonic trophoblast cells. hCG has several roles at the interface between the mother and the placenta. These roles include stimulating the ovarian corpus luteum to secrete progesterone, maintaining myometrial quiescence, and contributing to maternal immunotolerance towards the embryo. In addition, small amounts of hCG are produced in other tissues, including pituitary, testis, and colon. hCG is also expressed to some extent by several tumors, such as germ cell tumors and trophoblastic cancers, but its exact role in these cancers remains unknown.^{4–7}

hCG is a 37.5 kDa heterodimer consisting of two noncovalently linked subunits: hCG α , an alpha subunit identical in all glycoprotein hormones, and hCG β , a 23.5 kDa subunit unique to hCG that confers its biological specificity. Both hCG subunits are highly glycosylated; the carbohydrate chains attached to the possible glycosylation sites of the molecule make up about 30 percent of the hCG mass. The carbohydrate content of hCG varies substantially between different cell types, with advancing pregnancy, and with tumor progression. Due to this heterogeneity, hCG occurs in biological fluids in different forms, including intact hCG, free hCG β -subunit, free hCG α -subunit, and other nicked or truncated forms of the molecule. The biological effects of hCG are mediated through its binding to the receptor that recognizes hCG β , LHCGR, which is also activated by luteinizing hormone (LH). Both subunits are required for the biological activity of hCG.^{4,7–10}

The primary clinical use of hCG immunoassays is pregnancy detection. hCG can be detected from maternal blood already during the second week of pregnancy. Urine can also be used for hCG testing. Pregnancy tests are based on antibodies recognizing intact hCG, hCG β , or both. Additionally, measurement of hCG or its different molecular forms can be utilized for various other clinical purposes, such as pregnancy monitoring, diagnosis of pregnancy-related disorders, prenatal screening (e.g. for Down syndrome), as well as screening for trophoblastic tumors.^{47,8}

Medix Biochemica has more than three decades of experience in producing hCG mouse monoclonal antibodies for different diagnostic purposes. Currently, the company's products include eight anti-hCG antibodies recognizing either the hCG complex, free hCG β , or both. The antigen epitopes of hCG recognized by Medix Biochemica antibodies have been extensively characterized.⁸

Scientific publications: page 18

hCG antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
5004	100004	1	36	lgG ₁	ELISA, CLIA
5006	100005	5	24	lgG ₁	ELISA, CLIA
5008	100006	5	36	lgG ₁	ELISA, CLIA
5009	100008	5	18	lgG ₁	ELISA
5011	100009	1	24	lgG ₁	ELISA
5012	100368	5	36	lgG ₁	ELISA, CLIA
5014	100011	5	36	lgG ₁	ELISA
5016	100013	5	24	lgG ₁	ELISA, CLIA

Anti-human hCG monoclonal antibodies

Pair recommendations

		Detection										
					hCGI	beta				alpha subunit		
		5004	5006	5008	5009*	5011	5012**	5014	5016	5501	5503	6601
	5004	-	-	-	+	+	-	+	-	+	+ ^A	+ ^A
	5006	-	-	-	-	-	-	+	-	+	+	+
	5008	-	-	-	+	-	-	+ ^C	-	+	+	+
ture	5009*	+	+ ^A	+	-	-	-	+^	+	-	-	+
Capture	5011	+	+	+	-	-	-	+	+	-	-	+
	5012**	+	+ ^B	+	-	-	-	+	+	-	-	-
	5014	+	+	+	+	+	-	-	+	+	+	+
	5016	-	-	-	-	-	-	+	-	+	+	+

* Recognizes an epitope at the junction of alpha and beta subunits in hCG
** Recognizes only free beta subunit of hCG
A) Recommendation for detecting intact hCG molecule (alpha/beta pair)
B) Recommendation for detecting free beta subunit of hCG
C) Recommendation for detecting total hCG (both intact and free beta subunit of hCG)

Antibodies recognizing alpha subunit are recommended only for detection antibodies, since they also recognize other gonadotropins (including hCG, FSH, LH, and TSH).

hCG	Association rate constant,		Dissocia	ition rate constant,	Dissocia	Dissociation constant,		
antibody	k _{on} (1/Ms)		k _{off} (1/s)		K _D (M)	K _D (M)		
5004	hCG:	2.0 x 10 ⁶	hCG:	3.7 x 10 ⁻⁵	hCG:	1.9 x 10 ⁻¹¹ = 0.02 nM		
	hCGB:	5.0 x 10 ⁵	hCGB:	1.2 x 10 ⁻⁴	hCGB:	2.3 x 10 ⁻¹⁰ = 0.23 nM		
5006	hCG:	5.7 x 10 ⁶	hCG:	2.9 x 10 ⁻⁴	hCG:	5.0 x 10 ⁻¹¹ = 0.05 nM		
	hCGB:	2.3 x 10 ⁶	hCGB:	4.8 x 10 ⁻⁴	hCGB:	2.1 x 10 ⁻¹⁰ = 0.21 nM		
5008	hCG:	4.0 x 10 ⁶	hCG:	1.4 x 10 ⁻⁴	hCG:	3.4 x 10 ⁻¹¹ = 0.03 nM		
	hCGB:	1.6 x 10 ⁶	hCGB:	3.6 x 10 ⁻⁴	hCGB:	2.3 x 10 ⁻¹⁰ = 0.23 nM		
5009	hCG:	1.1 x 10⁵	hCG:	3.1 x 10 ⁻⁴	hCG:	2.8 x 10 ⁻¹⁰ = 0.28 nM		
	hCGB:	N/A	hCGB:	N/A	hCGB:	N/A		
5011	hCG:	2.4 x 10 ⁶	hCG:	4.4 x 10 ⁻⁴	hCG:	1.8 x 10 ⁻¹⁰ = 0.18 nM		
	hCGB:	1.9 x 10 ⁶	hCGB:	1.1 x 10 ⁻³	hCGB:	5.7 x 10 ⁻¹⁰ = 0.57 mM		
5012	hCG:	N/A	hCG:	N/A	hCG:	N/A		
	hCGB:	2.1 x 10 ⁴	hCGB:	1.7 x 10 ⁻⁴	hCGB:	8.1 x 10 ^{.10} = 0.81 nM		
5014	hCG:	2.5 x 10⁵	hCG:	1.5 x 10 ⁻⁴	hCG:	6.1 x 10 ⁻¹⁰ = 0.61 nM		
	hCGB:	2.8 x 10⁵	hCGB:	2.8 x 10 ⁻⁴	hCGB:	1.0 x 10 ⁻⁹ = 1.0 nM		
5016	hCG:	1.5 x 10 ⁶	hCG:	9.5 x 10 ⁻⁴	hCG:	6.4 x 10 ⁻¹⁰ = 0.64 nM		
	hCGB:	1.9 x 10 ⁶	hCGB:	5.2 x 10 ⁻³	hCGB:	2.8 x 10 ⁻⁹ = 2.8 nM		

Hormones

Epitopes and specificity

Berger et al. 2013⁸

Product	Epitope	hCG	hCGβ	hCGβcf	hCGn	hCGβn	-CTP hCG	-CTP- hCGβ	ISOBM Ab codes
5004	β2	✓	✓	✓	✓	✓	✓	✓	402
5006	β ₂	\checkmark	\checkmark	\checkmark	✓	✓	✓	✓	388
5008	β ₂	✓	\checkmark	✓	✓	✓	✓	✓	390
5016	β ₂	\checkmark	✓	✓	✓	✓	✓	✓	408
5012*	β ₇		\checkmark	\checkmark					409
5014	β,	✓	\checkmark		✓	✓	✓	✓	394
5011	β	\checkmark	✓		✓	✓	N/D	N/D	406
5009**	C ₂	\checkmark							387
5501	α	✓			✓		✓		383
5503	α	✓			✓		✓		385
6601	α	\checkmark			\checkmark		\checkmark		391

The darker turquoise color indicates strong reactivity.

Free beta subunit specific Mab


*) **) $\alpha\beta$ -heterodimer-specific Mab

N/D not determined

hCGn hCGbn -CTPhCG -CTPhCGβ

Nicked $\alpha\beta$ heterodimer, nicks in the region of aa hCG β 44-48 Nicked hCG β , nicks in the region of aa hCG β 44-48 hCG β truncated core hCG, missing most of the hCG β CTP (aa hCG β 121-145) hCG β truncated core hCG β (aa hCG β 1-120), missing most of the hCG β CTP

Epitope groups of hCG antibodies



hCG antibodies can be divided into different epitope groups. Antibodies in the same group detect the same or overlapping epitopes.

LH

Human luteinizing hormone (LH), or lutropin, is a glycoprotein hormone produced by the anterior pituitary gland. LH is produced in a pulsatile manner in response to gonadotropinreleasing hormone from the hypothalamus. LH is essential for normal reproduction in both males and females.¹¹⁻¹²

In females, LH together with follicle-stimulating hormone (FSH) regulates normal ovarian function. During the menstrual cycle, FSH initiates follicle growth leading to a rise in estrogen levels that triggers a peak in LH production, also known as the LH surge. The LH surge results in ovulation, and also stimulates corpus luteum formation and progesterone production, thus preparing the endometrium for possible pregnancy initiation. Furthermore, LH stimulates ovarian theca cells to produce androstenedione, which is converted to the female sex hormone estradiol. In males, LH stimulates testosterone production by the Leydig cells of the testis and is thus responsible for development of puberty, male secondary sexual characteristics, and spermatogenesis.^{11–15}

Like other glycoprotein hormones, LH is a heterodimer composed of a common 14 kDa alpha subunit, as well as a unique beta subunit that defines the hormone's functional specificity. The subunits of the 28 kDa LH are noncovalently linked to each other, and both subunits are highly glycosylated. Similar to hCG, the biological functions of LH are mediated through LHCGR binding. However, although LH and hCG activate the same receptor, their molecular properties and expression patterns differ significantly due to their different physiological roles. For example, the serum half-life of LH is less than 120 minutes, while it is several hours for hCG.^{13,16}

The clinical indications of assaying LH levels in blood or urine include predicting the timing of ovulation, investigating menstrual irregularities, and diagnosing early or delayed puberty. The LH test is often used in conjunction with other hormonal tests, such as tests for FSH or TSH levels.¹⁷⁻²⁰

Medix Biochemica has over three decades of experience in producing high-quality monoclonal anti-LH antibodies. Currently, the company offers four different mouse monoclonal antibodies for LH detection. Two have no crossreactivity with hCG, TSH, or FSH (clones 5303 and 5304), one has negligible cross-reactivity with TSH and FSH (clone 5301), and one has cross-reactivity to hCG (clone 5302). The association and dissociation constants of the antibodies have been measured using SPR.

Scientific publications: page 18

Anti-human LH monoclonal antibodies

LH antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
5301	100016	5	24	lgG ₁	ELISA
5302	100018	1	36	lgG ₁	ELISA
5303	100588	5	36	lgG ₁	ELISA
5304	100022	5	36	IgG ₁	ELISA

Hormones

Pair recommendations

		Detection							
			L	н		alpha subunit			
		5301	5302	5303	5304	5501			
	5301	-	+	+	+	+			
Capture	5302	+	-	-	-	+			
Capl	5303	+	-	-	-	-			
	5304	+	-	-	-	-			

LH antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _p (M)
5301	5.8 x 10 ⁶	1.3 x 10 ⁻⁵	2.3 x 10 ⁻¹² = 2.3 pM
5302	1.3 x 10 ⁶	2.6 x 10 ⁻⁵	2.1 X 10 ⁻¹¹ = 21 pM
5303	5.4 x 10 ⁶	3.4 x 10 ⁻⁵	6.3 X 10 ⁻¹² = 6.3 pM
5304	4.7 x 10 ⁶	1.6 x 10 ⁻⁵	3.4 x 10 ⁻¹² = 3.4 pM

FSH

Similar to LH, follicle-stimulating hormone (FSH) is a glycoprotein hormone secreted from the anterior pituitary gland under the control of gonadotropin-releasing hormone. FSH is an α/β heterodimer consisting of the glycoprotein hormone alpha subunit and a hormone-specific FSH β chain. FSH β is a 119 amino acid peptide with two N-linked glycosylation sites. FSH mediates its biological functions by binding to the FSH receptor (FSHR), which is expressed mainly by ovarian granulosa cells and testicular Sertoli cells, but also by osteoclasts and certain tumor cells.^{11,21–23}

FSH has several essential functions related to mammalian development and reproduction. In both sexes, FSH stimulates germ cell maturation. During male development, FSH promotes spermatogonia mitosis and ensures normal testes development by activating proliferation of Sertoli cells. In adult males, FSH maintains normal spermatogenesis and Sertoli cell activities important for testis function and maintenance of male fertility. In females, FSH drives follicle growth and granulosa cell estrogen production. FSH also prevents granulosa cell apoptosis and stimulates their proliferation. Furthermore, there is a small rise in FSH levels at the end of the menstrual cycle luteal phase; this seems to be relevant for the initiation of the next ovulatory cycle. Besides its reproductive roles, FSH has also been associated with postmenopausal bone loss and neovascularization surrounding malignant tumors; however, these functions require further investigation.^{11,15,22,24}


Clinically, FSH is commonly used for ovarian stimulation in women who are undergoing infertility treatment or oocyte retrieval. FSH can also be used for treating anovulatory infertility in women or hypogonadotropic hypogonadism in men. LH is often used as a supportive treatment with FSH therapy. FSH tests are often carried out together with an LH test. Assaying serum or urine FSH levels can be used to investigate female infertility, menstrual irregularities, or to diagnose conditions related to ovarian or testicular dysfunction.^{15,22,25}

Medix Biochemica offers one FSH-specific mouse monoclonal antibody (6602) that can be used as a pair with an alpha subunit antibody (6601) in diagnostic tests.

Anti-human FSH monoclonal antibody

FSH	Product	Concentration	Shelf life	Subclass	Applications
antibody	code	(mg/mL)	(months at +2–8°C)		tested
6602	100067	5	24	lgG ₁	ELISA

Pair recommendations



* alpha subunit

FSH	Dissociation constant,
antibody	K _D (M)
6602	3.0 x 10 ⁻¹¹ = 0.03 nM

Alpha subunit

The alpha subunit is common to all glycoprotein hormones, also known as gonadotropins. The glycoprotein hormones are specific to vertebrates and include hCG, FSH, LH, and TSH. In humans, the alpha subunit consists of 92 amino acids and has an average molecular weight of 14 kDa. In all gonadotropins, the alpha subunit binds non-covalently to the beta subunit forming a heterodimer; each beta subunit is hormone-specific and mediates receptor binding and the specific functions of the glycoprotein hormones. However, the alpha subunit is also required for the biological activity of these hormones.^{5,10,26} Medix Biochemica offers three mouse monoclonal antibodies (5501, 5503, and 6601) for the detection of the human glycoprotein hormone alpha subunit. The association and dissociation characteristics of the binding reaction between the alpha subunit antibodies and the alpha subunit have been measured using multiplexed SPR technology. Anti-human alpha subunit specific antibodies are recommended only as detection antibodies with human glycoprotein hormone (hCG, LH, FSH and TSH) beta subunit specific antibodies. Antihuman alpha subunit antibodies bind to several glycoprotein hormones and therefore are not recommended to be used as capture antibodies.


Anti-human alpha subunit monoclonal antibodies

Alpha subunit antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
5501	100035	1	36	lgG ₁	ELISA
5503	100037	1	6	lgG ₁	ELISA
6601	100066	5	36	lgG _{2a}	ELISA

Kinetic parameters

Alpha subunit antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, $K_{\rm p}$ (M)
5501	1.9 x 10 ⁶	1.1 x 10 ⁻³	5.7 x 10 ⁻¹⁰ = 0.57 nM
5503	1.7 x 10 ⁶	5.3 x 10 ⁻⁴	3.2 x 10 ⁻¹⁰ = 0.32 nM
6601	8.9 x 10 ⁵	3.8 x 10 ⁻⁵	4.2 x 10 ⁻¹¹ = 0.04 nM

Binding of glycoprotein hormones



Glycoprotein hormones (TSH, FSH, LH, and hCG) have an identical alpha subunit (α) and a hormonespecific beta subunit (β).

SHBG

Sex hormone binding globulin (SHBG) is a steroid-binding serum glycoprotein synthesized in the liver. Synthesized SHBG is secreted into the circulation, where it binds the sex hormones testosterone and estradiol with high affinity and specificity. Thus, SHBG regulates sex hormone bioavailability and affects steroid responses by inhibiting their function. Some SHBG is also expressed by testicular germ cells, suggesting that the protein may have influence on sperm function. SHBG serum levels are regulated primarily by sex hormone blood levels, but also e.g. by nutritional factors.^{27,28}

SHBG is an approximately 90 kDa homodimeric glycoprotein that is present in blood only as dimer. Dimerization occurs automatically upon synthesis and secretion of the molecule, and is promoted by high-affinity steroid ligands of the hormone. Upon dimerization, eight hydrogen bonds form between β -sheets of the monomers. This leads to an elongated, cylinder-shaped mature dimer. ^27,28

SHBG immunoassays are often used to assess patient testosterone status, such as in male infertility treatment or female hirsutism. In addition, several conditions affect SHBG plasma levels either by increasing or decreasing it. For example, low SHBG levels can indicate diabetes, hypothyroidism, or polycystic ovary syndrome. High SHBG levels may be associated with pregnancy, hyperthyroidism, or anorexia nervosa. Furthermore, recent research has also linked high SHBG levels with breast and testicular cancer. Epidemiological studies indicate that low serum SHBG levels predict a higher risk for metabolic syndrome.²⁷²⁹

Medix Biochemica's selection of mouse monoclonal anti-SHBG antibodies includes three products (6001, 6002, and 6007).

Anti-human SHBG monoclonal antibodies

SHBG antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
6001	100539	5	24	lgG ₁	ELISA
6002	100540	5	24	lgG ₁	ELISA
6007	100569	5	24	lgG ₁	ELISA

Pair recommendations

			Detection	1
		6001	6002	6007
ture	6001	-	+	+
Þt	6002	-	-	-
ů	6007	+	-	-

SHBG antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K_{D} (M)
6001	4.1 x 10 ⁵	4.7 x 10 ⁻⁵	1.2 x 10 ⁻¹⁰ = 0.12 nM
6002	1.5 x 10⁵	1.0 x 10 ⁻⁴	6.8 x 10 ⁻¹⁰ = 0.68 M
6007	2.9 x 10⁵	1.1 × 10 ⁻⁴	3.7 x 10 ⁻¹⁰ = 0.37 nM

Progesterone

Progesterone (PR) is a steroid hormone produced in the corpus luteum after ovulation and in the placenta during pregnancy. Progesterone is also produced in the adrenal glands. The main functions of progesterone include induction of endometrial transition from a proliferative to a secretory stage during the menstrual cycle, facilitating blastocyst nesting, as well as maintenance of pregnancy. Progesterone also promotes insulin release and may have some neuroprotective effects, including slowing down Alzheimer's disease progression.³⁰⁻³²

In women, progesterone levels are generally low (approximately <1 ng/mL) before ovulation and rise thereafter up to 10 to 35 ng/mL. During pregnancy, progesterone levels may reach 100 to 300 ng/mL. In early pregnancy, progesterone measurements are used in addition to hCG measurements for more precise diagnosis of an ectopic or otherwise dysfunctional pregnancy. As low progesterone levels can lead to pregnancy termination, progesterone levels are often monitored during high-risk pregnancies to evaluate and ensure fetal health.³⁰⁻³¹

Medix Biochemica offers four anti-progesterone mouse monoclonal antibodies (1801, 1802, 1803, and 1805). Due to the small size of the analyte, competitive immunoassay is recommended for progesterone level determination. Progesterone conjugated to BSA from C6 has been used as the immunogen for the production of Medix Biochemica antiprogesterone antibodies.

Anti-human progesterone monoclonal antibodies

PR antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
1801	100248	5	24	lgG ₁	ELISA
1802	100249	5	36	lgG _{2a}	ELISA
1803	100250	5	24	lgG _{2a}	ELISA
1805	100252	5	36	lgG _{2a}	ELISA

Prolactin

Prolactin (PRL) is a peptide hormone secreted primarily by the lactotroph cells of the pituitary gland. PRL is also produced by various extrapituitary cells and sites, including immune cells, neurons, prostate, mammary epithelium, and skin. During pregnancy, PRL acts in concert with other hormones, including progesterone and insulin, to promote mammary gland growth and development. In the postpartum period, PRL stimulates lactation by initiating milk synthesis and maintaining its secretion. PRL production reduces the secretion of gonadotropin-releasing hormone, thereby suppressing LH and FSH secretion and inhibiting ovulation. Besides its effects on pituitary-gonadal processes, PRL is involved in various other processes. These include regulating osmotic balance, increasing bile and insulin secretion, and regulating the immune system, for example via stimulation of T and NK cell cytokine expression.^{33,34}

PRL is a 23 kDa hormone that consists of 199 amino acids. Structurally, PRL resembles members of cytokine family molecules and consists of a single polypeptide chain with three intramolecular disulfide bonds. The various effects of PRL are mediated by its binding to a specific PRL receptor (PRLR), a type I cytokine receptor family member expressed in various tissues and as several different isoforms. PRL secretion is controlled mainly by dopamine, which inhibits the release of PRL. There are also several other molecules that influence PRL secretion either by stimulation (e.g. thyrotropin-releasing hormone, serotonin, oxytocin) or by inhibition (e.g. serotonin, histamine, noradrenaline).^{34,35}


PRL blood levels are generally relatively low, but high PRL levels are expected during pregnancy and nursing. High PRL levels in other situations may indicate a potential disorder. For example, high PRL levels may be related to pituitary or other PRL-releasing tumors, certain kidney or liver diseases, or hypothalamus diseases. PRL levels are often measured in cases of infertility or irregular menstrual periods, testicular dysfunction, or symptoms suggesting prolactinoma, including visual impairment. PRL levels vary significantly during the day, which should be taken into account when using them as a diagnostic tool.^{34,36-38}

Medix Biochemica offers two anti-PRL mouse monoclonal antibodies (5601 and 5602), which can be used as a pair in diagnostic tests.

Anti-human prolactin monoclonal antibodies

PRL antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
5601	100038	1	36	lgG ₁	ELISA
5602	100360	5	24	IgG ₁	ELISA

Pair recommendations



PRL antibody	Dissociation constant, K _D (M)
5601	3.8 x 10 ⁻¹¹ = 0.04 nM
5602	1.05 x 10 ⁻¹⁰ = 0.1 nM

T4

Thyroxine, or 3,5,3',5'-tetraiodothyronine (T4) is a small hapten prohormone synthesized and secreted by the thyroid gland. T4 is secreted in the bloodstream, where almost all of it is bound to three major serum transport proteins: mostly thyroxine-binding globulin (TBG), but also transthyretin and human serum albumin. Only 0.02 to 0.03% of total serum T4 circulates in the unbound form. The thyroid-hormone-binding plasma proteins are evenly distributed to maintain a stable T4 concentration in the body. The thyroid gland secretes an estimated 110 nmol of T4 daily. T4 synthesis is controlled by thyroid-stimulating hormone (TSH). Secreted T4 binds to the nuclear thyroid hormone receptor (TR), which is involved in the regulation of all intermediary metabolic processes of carbohydrates, lipids, and proteins. T4 also serves as a reservoir for the more metabolically active thyroid hormone, T3, which can be produced from T4 by deiodination. In the fetus and during childhood, thyroid hormones are critical for functions such as brain development, neuronal differentiation, and formation of neural processes.³⁹⁻⁴²

Medix Biochemica offers two anti-T4 mouse monoclonal antibodies (6901 and 6902), which can be used for a competitive immunoassay for T4 detection.

Anti-human T4 monoclonal antibodies

T4 antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
6901	100348	5	36	lgG ₁	ELISA
6902	100801	5	N/D	lgG ₁	ELISA

T4 antibody	Dissociation constant, K_{D} (M)
6901	1.0 x 10 ⁻¹⁰ = 0.1 nM
6902	N/D

Thyroglobulin

Thyroglobulin is the predominant protein produced by the thyroid gland. This homodimeric 660 kDa glycoprotein protein acts as a substrate for the synthesis of the thyroid hormones triiodothyroxine (T3) and thyroxine (T4). In addition, this T3- and T4-precursor serves as a storage protein for iodide that ensures sufficiency of this essential element.⁴³

Thyroglobulin is also produced by thyroid cancer cells, and is thus commonly used to monitor the treatment of cancer patients who have undergone thyroidectomy. However, as thyroglobulin is also produced by normal thyroid cells, it is not suitable as a diagnostic cancer biomarker. Blood thyroglobulin levels can also be elevated in some cases of Grave's disease. Thyroglobulin has recently been suggested a potential biomarker for iodine deficiency.⁴²⁻⁴⁵

Medix Biochemica's antibody selection contains four antithyroglobulin mouse monoclonal antibodies (2802, 2803, 2804, and 2805).

Anti-human thyroglobulin monoclonal antibodies

Thyroglobulin antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
2802	100331	5	18	lgG ₁	ELISA
2803	100332	5	36	IgG _{2b}	ELISA
2804	100333	5	36	IgG ₁	ELISA
2805	100334	5	36	lgG ₁	ELISA

Pair recommendations

		Detection			
		2802	2803	2804	2805
cure	2802	-	-	-	-
	2803	-	-	+	+
Captur	2804	-	-	-	-
-	2805	+	+	+	-

Thyroglobulin antibody	Association rate constant, k _{on} (1/Ms)	Dissociation rate constant, k _{off} (1/s)	Dissociation constant, K _D (M)
2802	8.6 x 10 ³	1.2 x 10 ⁻⁵	1.4 x 10 ⁻⁹ = 1.4 nM
2803	1.7 x 10⁵	3.7 x 10 ⁻⁵	2.2 x 10 ⁻¹⁰ = 0.22 nM
2804	7.2 x 10 ⁴	7.9 x 10 ⁻⁴	1.1 x 10 ⁻⁸ = 11 nM
2805	6.2 x 10 ⁵	1.3 x 10 ⁻⁵	2.1 x 10 ⁻¹¹ = 0.02 nM

Growth hormone


Human growth hormone (GH) is a single-chain 22 kDa peptide hormone synthesized and secreted by the anterior pituitary gland. In accordance with its name, GH stimulates various bodily processes related mostly to growth and regeneration. Besides stimulating longitudinal bone growth, GH regulates cell growth and differentiation, as well as metabolic processes in various tissues, such as liver, adipose tissue, muscle, and bone. GH release is regulated primarily by growth hormone releasing hormone (GHRH), which induces GH secretion, and by somatostatin, which inhibits it. While males secrete GH in a pulsatile manner, females have been reported to exhibit more continuous GH secretion. GH function is mediated by its binding to the GH receptor, a type I cytokine receptor expressed primarily in the liver. Clinically, GH testing is mainly used for identifying growth hormone deficiency or, less commonly, excess GH related to acromegaly and gigantism.⁴⁶⁻⁴⁹

Medix Biochemica's selection of mouse monoclonal antibodies includes two anti-GH antibodies (5801 and 5802). These can be used as a pair in diagnostic tests.

Anti-human growth hormone monoclonal antibodies

GH antibody	Product code	Concentration (mg/mL)	Shelf life (months at +2–8°C)	Subclass	Applications tested
5801	100041	1	36	lgG ₁	ELISA
5802	100042	1	36	lgG ₁	ELISA

Pair recommendations



GH antibody	Dissociation constant, K _p (M)
5801	5.0 x 10 ⁻¹¹ = 0.05 nM
5802	1.0 x 10 ⁻¹⁰ = 0.1 nM

Page 3

Scientific articles describing the use of Medix Biochemica's monoclonal anti-TSH antibodies include:

- von Lode P, Hagren V, Palenius T & Lövgren T (2003). One-step quantitative thyrotropin assay for the detection of thyrotropin assay for the detection of hypothyroidism in point-of-care conditions. Clin Biochem 36:121–128.
- Wu FB, Han SQ & He YF (2002). Time-resolved immunofluorometry of serum hTSH with enhanced sensitivity. J Immunoassay Immunochem 23:191–210.
- Helenius T & Tikanoja S (1986). A sensitive and practical immunoradiometric assay of thyrotropin. Clin Chem 32:514–518.

Page 5

Scientific articles describing the use of Medix Biochemica's monoclonal anti-hCG antibodies include:

- Berger P, Paus E, Hemken PM et al. (2013). Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7 workshop. Tumour Biol 34:4033–4057.
- von Lode P, Rainaho J & Pettersson K (2004). Quantitative, wide-range, 5-minute point-of-care immunoassay for total human chorionic gonadotropin in whole blood. Clin Chem 50:1026–1035.
- Stenman UH, Tanner P, Ranta T et al. (1982). Monoclonal antibodies to chorionic gonadotropin: use in a rapid radioimmunoassay for gynecologic emergencies. Obstet Gynecol 59:375–377.

Page 8

Publications describing the use of Medix Biochemica's monoclonal LH antibodies:

- Bielmeier SR, Best DS & Narotsky MG (2004). Serum hormone characterization and exogeneous hormone rescue of bromodichloromethane-induced pregnancy loss in the F344 rat. Toxicol Sci 77:101–108.
- Manna PR, Joshi L, Reinhold VN et al. (2002). Synthesis, purification and structural and functional characterization of recombinant form of a common genetic variant of human luteinizing hormone. Hum Mol Gen 11:301–315.
- Yano K, Hashida S, Okuno A & Ishikawa E (1988). Serum concentration of chorionic gonadotropin-like substance and luteinizing hormone in children measured by sensitive enzyme immunoassay. J Clin Lab Anal 2:220–224.
- Federici MM, Fraser R, Lundqvist C & Lankford JC (1982). Production and characterization of monoclonal antibodies human luteinizing hormones. Fed Proc 41:596, Abstract 1907.

References:

- Grossmann M, Szkudlinski MW, Wong R et al. (1997). Substitution of the seat-belt region of the thyroid-stimulating hormone (TSH) beta-subunit with the corresponding regions of choriogonadotropin or follitropin confers luteotropic but not follitropic activity to chimeric TSH. J Biol Chem 272:15532– 15540.
- Szkudlinski MW, Fremont V, Ronin C & Weintraub BD (2002). Thyroid-stimulating hormone and thyroid-stimulating hormone receptor structure-function relationships. Physiol Rev 82:473–502.
- Thienpont LM, van Uytfanghe K, van Houcke S et al. (2014). A progress report of the IFCC committee for standardization of thyroid function tests. Eur Thyroid J 3:109–116.
- 4. Cole LA (2010). Biological functions of hCG and hCG-related molecules. Reprod Biol Endocrinol 8:102.
- Fournier T, Guibourdenche J & Evain-Brion D (2015). Review: hCGs: different sources of production, different glycoforms and functions. Placenta 36 Suppl 1:S60–5.
- Montagnana M, Trenti T, Aloe R et al. (2011). Human chorionic gonadotropin in pregnancy diagnostics. Clin Chim Acta 412:1515–1520.
- Stenman UH & Alfthan H (2013). Determination of human chorionic gonadotropin. Best Pract Res Clin Endocrinol Metab 27:783–793.
- Berger P, Paus E, Hemken PM et al. (2013). Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7 workshop. Tumour Biol 34:4033–4057.
- Lapthorn AJ, Harris DC, Littlejohn A et al. (1994). Crystal structure of human chorionic gonadotropin. Nature 369:455– 461.
- Stenman UH, Tiitinen A, Alfthan H & Valmu L (2006). The classification, functions and clinical use of different isoforms of HCG. Hum Reprod Update 12:769–784.
- 11. Richards JS & Pangas SA (2010). The ovary: basic biology and clinical implications. J Clin Invest 120:963–972.
- 12. Narayan P (2015). Genetic models for the study of luteinizing hormone receptor function. Front Endocrinol (Lausanne) 6:152.
- Nilsson C, Seppälä M & Pettersson K (2001). Immunological characterization of human luteinizing hormone with special regard to a common genetic variant. J Endocrinol 168:107–116.
- Raju GA, Chavan R, Deenadayal M et al. (2013). Luteinizing hormone and follicle stimulating hormone synergy: A review of role in controlled ovarian hyper-stimulation. J Hum Reprod Sci 6:227–234.
- Reed BG & Carr BR (2015). The normal menstrual cycle and the control of ovulation. In: Endocrinology of Female Reproduction (ed. Rebar R). Available in: http://www.endotext.org/chapter/ the-normal-menstrual-cycle-and-the-control-of-ovulation/ (Accessed 05/2016).

- 16. Casarini L, Lispi M, Longobardi S et al. (2012). LH and hCG action on the same receptor results in quantitatively and qualitatively different intracellular signalling. PLoS One 7:e46682.
- Zollner U (2011). Hormonal diagnosis of menstrual irregularities or secondary amenorrhoea. In: Amenorrhea. Available in: http:// www.intechopen.com/books/amenorrhea (Accessed 06/2016).
- Corsan GH, Ghazi D & Kemmann E (1990). Home urinary luteinizing hormone immunoassays: clinical applications. Fertil Steril 53:591–601.
- Houk CP, Kunselman AR & Lee PA (2009). Adequacy of a single unstimulated luteinizing hormone level to diagnose central precocious puberty in girls. Pediatrics 123:e1059–63.
- Chin VL, Cai Z, Lam L et al. (2015). Evaluation of puberty by verifying spontaneous and stimulated gonadotropin values in girls. J Pediatr Endocrinol Metab 28:387–392.
- Fan QR & Hendrickson WA (2005). Structure of human folliclestimulating hormone in complex with its receptor. Nature 433:269–277.
- Laan M, Grigorova M & Huhtaniemi IT (2012). Pharmacogenetics of follicle-stimulating hormone action. Curr Opin Endocrinol Diabetes Obes 19:220–227.
- Jiang X, Liu H, Chen X et al. (2012). Structure of folliclestimulating hormone in complex with the entire ectodomain of its receptor. Proc Natl Acad Sci U S A 109:12491–12496.
- Busch AS, Hagen CP, Almstrup K et al. (2016). Genetic variations altering FSH action affect circulating hormone levels as well as follicle growth in healthy peripubertal girls. Hum Reprod 31:897–904.
- Simoni M & Casarini L (2014). Mechanisms in endocrinology: Genetics of FSH action: a 2014-and-beyond view. Eur J Endocrinol 170:R91–107.
- Fiddes JC & Goodman HM (1981). The gene encoding the common alpha subunit of the four human glycoprotein hormones. J Mol Appl Genet 1:3–18.
- 27. Simo R, Saez-Lopez C, Barbosa-Desongles A et al. (2015). Novel insights in SHBG regulation and clinical implications. Trends Endocrinol Metab 26:376–383.
- Thaler MA, Seifert-Klauss V & Luppa PB (2015). The biomarker sex hormone-binding globulin – from established applications to emerging trends in clinical medicine. Best Pract Res Clin Endocrinol Metab 29:749–760.
- Sheehan MT (2004). Polycystic ovarian syndrome: diagnosis and management. Clin Med Res 2:13–27.
- Taraborrelli S (2015). Physiology, production and action of progesterone. Acta Obstet Gynecol Scand 94 Suppl 161:8–16.
- Lessey BA, Killam AP, Metzger DA et al. (1988). Immunohistochemical analysis of human uterine estrogen and progesterone receptors throughout the menstrual cycle. J Clin Endocrinol Metab 67:334–340.

- Wang JM, Irwin RW, Liu L et al. (2007). Regeneration in a degenerating brain: potential of allopregnanolone as a neuroregenerative agent. Curr Alzheimer Res 4:510–517.
- Shelly S, Boaz M & Orbach H (2012). Prolactin and autoimmunity. Autoimmun Rev 11:A465–70.
- Ignacak A, Kasztelnik M, Sliwa T et al. (2012). Prolactin--not only lactotrophin. A "new" view of the "old" hormone. J Physiol Pharmacol 63:435–443.
- Brooks CL (2012). Molecular mechanisms of prolactin and its receptor. Endocr Rev 33:504–525.
- 36. Ress C, Maeser PA, Tschoner A et al. (2014). Serum prolactin in advanced chronic liver disease. Horm Metab Res 46:800–803.
- Serri O, Chik CL, Ur E & Ezzat S (2003). Diagnosis and management of hyperprolactinemia. CMAJ 169:575–581.
- Jenkins PJ & Besser M (2003). Disorders of the anterior pituitary. In: Oxford Textbook of Medicine (ed. Warrell D, Cox T, Firth JD & Benz EJ): 191–203. Oxford Medical Publications, Oxford.
- Hulbert AJ (2000). Thyroid hormones and their effects: a new perspective. Biol Rev Camb Philos Soc 75:519–631.
- 40. Moreno M, de Lange P, Lombardi A et al. (2008). Metabolic effects of thyroid hormone derivatives. Thyroid 18:239–253.
- Pharoah P, Buttfield IH & Hetzel BS (2012). Neurological damage to the fetus resulting from severe iodine deficiency during pregnancy. Int J Epidemiol 41:589–592.
- American Thyroid Association (2014). Thyroid Function Tests. Available in: http://www.thyroid.org/wp-content/uploads/ patients/brochures/FunctionTests_brochure.pdf (Accessed 05/2016).
- 43. Di Jeso B & Arvan P (2016). Thyroglobulin from molecular and cellular biology to clinical endocrinology. Endocr Rev 37:2–36.
- Ma ZF & Skeaff SA (2014). Thyroglobulin as a biomarker of iodine deficiency: a review. Thyroid 24:1195–1209.
- 45. van de Graaf SA, Ris-Stalpers C, Pauws E et al. (2001). Up to date with human thyroglobulin. J Endocrinol 170:307–321.
- Chia DJ (2014). Minireview: mechanisms of growth hormonemediated gene regulation. Mol Endocrinol 28:1012–1025.
- Vijayakumar A, Yakar S & Leroith D (2011). The intricate role of growth hormone in metabolism. Front Endocrinol (Lausanne) 2:32.
- Melmed S, Casanueva FF, Klibanski A et al. (2013). A consensus on the diagnosis and treatment of acromegaly complications. Pituitary 16:294–302.
- 49. Davies JH & Cheetham T (2014). Investigation and management of tall stature. Arch Dis Child 99:772–777.

More information on our products **medixbiochemica.com**



Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com

CPS = Counts per second CLIA = Chemiluminescence immunoassay ELISA = Enzyme-linked immunosorbent assay FIA = Fluoroimmunoassay IT = Immunoturbidimetry LF = Lateral flow N/A = Not Applicable N/D = Not Determined

Copyright © 6/2020 Medix Biochemica. All rights reserved.

MedixMAB and MedixAntigens are trademarks of Medix Biochemica. Medix Biochemica reserves the right to make changes and improvements to any of the products described in this document without prior notice.