




# Application data (LipiRADICAL)


## Cell-based imaging

Hepa1-6 cells were treated with "LipiRADICAL Green". For inducing an LPO signal, the cells were co-treated with diethylnitrosamine (DEN) and "LipiRADICAL Green", an LPO initiator. Immediately after DEN addition, the cells were observed by confocal microscopy with 2 min interval. The fluorescent signal of "LipiRADICAL Green" from the DEN-treated cells clearly increased.



## in vitro detection of lipid radicals derived from LDL

Purified low-density lipoprotein (LDL) was mixed with pro-oxidants hemin or AAPH and "LipiRADICAL Green" and the green fluorescence was measured. Both hemin and AAPH increased green fluorescence indicating the production of lipid radicals from LDL particles in a time-dependent manner.





## Application data (OH-Pen)

### Inhibition of nitrosamine-induced carcinogenesis by OH-Pen

Rats received diethylnitrosamine (DEN, 100 mg/kg body weight), which is a well-known hepatic procarcinogen. Subsequently, rats received OH-Pen (2.5 μmol/kg body weight) by intraperitoneal injection after 1 hour DEN administration. For the acute model and chronic model, livers were dissected after 24 hours and 12 weeks DEN administration, respectively. In all panels, OH-Pen clearly suppressed DEN-induced hepatocellular carcinoma.

- A. Livers from chronic hepatocellular carcinoma model and total foci number.
- B. Quantification of LPO-derived aldehydes in acute model livers.
- C. Quantification of tissue damage markers.



### What is Lipid Peroxidation (LPO)?

*Memo*

Lipid peroxidation (LPO) is one of the several degradation processes of lipids under oxidative stress. In the termination reaction, antioxidants donate a hydrogen atom to the lipid peroxy radical (LOO·) species resulting in the formation of many different aldehydes including malondialdehyde (MDA), acrolein, propanal, hexanal, and 4-hydroxynonenal (4-HNE). These reactive aldehydes are considered as causative factors of organ injury, ferroptosis and ER stress.

<Manufacturer: FNA>

Product Name	Code	Size	Price
LipiRADICAL Green <Lipid Radical Detection Reagent>	FDV-0042	0.1 mg	

<Manufacturer: FNA>

Product Name	Code	Size	Price
OH-Pen <Lipid Radical Inhibitor>	FDV-0043	0.1 mg	

※ All products here are research use only, not for diagnostic use.  
 ※ Specs might be changed for improvement without notice.  
 ※ Numbers after "#" represents product code.

※ Company name and product name are trademark or registered mark.  
 ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor

2BScientific Ltd,



Cherwell Innovation Centre, 77 Heyford Park,  
 Upper Heyford, OX25 5HD, UK  
 Phone: +44 (0)1869 238033  
 Fax : +44 (0)1869 238034  
 General: info@2BScientific.com

Funakoshi Co., Ltd. (JAPAN)

Phone: +81-3-5684-6296  
 Mail: export@funakoshi.co.jp



<Contact/



/Newsletter Signup/>