

JAPAN MADE QUALITY

## DNs-Rh


For more information : [https://www.funakoshi.co.jp/exports\\_contents/81304](https://www.funakoshi.co.jp/exports_contents/81304)

**DNs-Rh** is a novel fluorogenic substrate for GSTs (Glutathione S-Transferases). Compared to conventional reagents, DNs-Rh can be used for live cell imaging with high specificity and sensitivity.

This product has been commercialized under the license from Nagoya University.

| Probe            |                   | DNs-Rh             | Conventional method (CDNB) |
|------------------|-------------------|--------------------|----------------------------|
| Detection Method |                   | Green Fluorescence | UV                         |
| Application      | Live Cell Imaging | Yes                | No                         |
|                  | Flow Cytometry    | Yes                | No                         |
| Sensitivity      |                   | High               | Low                        |
| Throughput       |                   | High               | Low                        |
| Specificity      |                   | High               | Low                        |
| Co-staining      |                   | Yes                | No                         |

### Principle



### Application


- Intracellular GST activity assay in live cells (Imaging or flow cytometry)  
**Note : Commercial FITC filter sets are compatible.**
- *in vitro* GST activity assay  
**Note : Broad specificity for GST family**

### Product Information

[ Manufacturer : FNA ]

| Product Name                                   | Code     | Size     | Storage |
|--|----------|----------|---------|
| DNs-Rh <Cell-based GST Activity Assay Reagent> | FDV-0030 | 0.1 μmol | -20°C   |

## Example Data





### ▲ Fig.1 *in vitro* GST activity assay

Assay solution containing 1  $\mu$ M DNs-Rh, 1 mM GSH and 10  $\mu$ g/ml recombinant human GSTP1-1 was incubated at 37°C for 30 min. DNs-Rh only showed no fluorescence. In the presence of both GSTP1-1 and GSH, fluorescence was dramatically increased and about 60% of DNs-Rh was converted to rhodamine 110 at 30 min. However, in the absence of GSTP1, fluorescent change by the free thiol of GSH was much slow (3% at 30 min).

### Fig.2 Monitoring intracellular GST activity in live cells ►

HeLa cells were treated with 2.5  $\mu$ M of DNs-Rh for 30 min. In the case of addition of CNBSF, a potent irreversible GST inhibitor, cells were pre-treated with 1 mM CNBSF for 15 min prior to DNs-Rh addition. After DNs-Rh incubation, the cells were washed and observed fluorescent microscopy (Ex. 480 nm / Em. 535 nm). CNBSF is also available from Funakoshi, catalog #FDV-0031 (see right).

NOTE: Using this probe, you can observe the intracellular localization of rhodamine 110 converted from DNs-Rh by GSTs. Fluorescent intensity is corresponding to cellular GST activity. However, localization is not equal to localization of GSTs.



[ Manufacturer : FNA ]

| Product Name                          | Code     | Size  | Storage |
|---------------------------------------|----------|-------|---------|
| CNBSF<br><Irreversible GST Inhibitor> | FDV-0031 | 10 mg | -20°C   |

## Background about GSTs and advantage of DNs-Rh

As many studies suggested, expression level of GSTs is significantly increased in cancer cells, GSTs are considered as anti-cancer drug-resistant enzymes in malignant cancer cells through the neutralization of drugs. To understand biological functions of GSTs, tools for monitoring GST activity are very important. Although several reagents including CDNB (conventional GSTs probe) for this purpose have been developed, no tool to measure intracellular GST activity was commercially available.

### DNs-Rh is a novel fluorogenic substrate of GSTs discovered by Dr. Hiroshi Abe, Nagoya University.

This probe emits very low fluorescence (quantum yield = 0.0007) at normal state. After de-protected by thiols, rhodamine 110 is released and exhibits strong fluorescence (quantum yield = 0.645, S/N ratio ~900). DNs-Rh can be used as a fluorogenic GST activity assay probe. An important advantage of this probe is high cell-permeability and this probe is applied to intracellular GST activity under live cell condition. As DNs group is a well characterized substrate for various types of GST members, DNs-Rh is able to monitor pan-GST activity both in cell and *in vitro*.

DNs-Rh is a powerful tool not only to investigate GST activity in live cell upon any biological stimulation, but also to develop GST inhibitors under live cell condition.

#### NOTE

※ All products here are research use only, not for diagnostic use.  
 ※ Specs might be changed for improvement without notice.

※ Company name and product name are trademark or registered mark.  
 ※ Please contact your local distributors for orders, quote request and inquiry.

Your Local Distributor

2BScientific Ltd,



Cherwell Innovation Centre, 77 Heyford Park,  
 Upper Heyford, OX25 5HD, UK  
 Phone: +44 (0)1869 238033  
 Fax : +44 (0)1869 238034  
 General: info@2BScientific.com

Funakoshi Co., Ltd.

Address: 9-7 Hongo 2-Chome, Bunkyo-ku,  
 Tokyo 113-0033 JAPAN  
 Phone : +81-3-5684-6296  
 Fax : +81-3-5684-6297  
 Email : export@funakoshi.co.jp

FNA-1905-F13 (2019.05)