



## FIBROBLASTS ISOLATION FROM HUMAN SKIN

### Enzymes

COL G recombinant collagenase class I ACTIVITY  $\geq 3.0$  Units/mg (Pz Grassmann)

COL H recombinant collagenase class II ACTIVITY  $\geq 30.0$  Units/mg (Pz Grassmann)

### Preparation of COL G and COL H stock solutions

1. Solution A: Reconstitute COL G in H<sub>2</sub>O, do not exceed 30 U/ml final concentration.
2. Solution B: Reconstitute COL H in H<sub>2</sub>O, do not exceed 300 U/ml final concentration.
3. Filter for sterility (0.22  $\mu$ m) and make aliquots of Solution A and Solution B.
4. Annotate the Unit/ml value of the stock solutions and store at -20°C

### Procedure

1. Transfer the tissue in a petri dish and wash with PBS (Dulbecco's Phosphate Buffered Saline Solution without Ca<sup>2+</sup> Mg<sup>2+</sup>). Cut the skin in 0.5 x 1.0 cm<sup>2</sup> fragments.
2. Add 5 ml of Dispase II and incubate for 3h at 37°C or overnight at 4°C.
3. Separate dermis from epidermis with sterile tweezers and transfer the dermis in another petri dish. Wash with PBS.
4. Mince the dermis and incubate the pieces in a petri dish with 15 ml of collagenases solution for about 4 h (checking regularly the dissociation proceeding) at 37°C containing 0,5 unit/ml of G and 5 unit/ml of H plus neutral protease (1 mg/ml) in PBS (Dulbecco's Phosphate Buffered Saline Solution with Ca<sup>2+</sup> Mg<sup>2+</sup>).
5. Collect digested dermis from the petri dish and add an amount of culture medium equal to the double of the collected volume.
6. Centrifuge 300x g for 4 min. Remove the supernatant and resuspend pellet in culture medium.
7. Seed a final volume of 5 ml in T25 flask and incubate at 37°C, 5% CO<sub>2</sub>.

**Note:** *This protocol is meant to be a starting point, all isolation procedures require an individual optimization. COL G: COL H concentration, protease addition and digestion time can be adjusted experimentally.*

If you have any question please contact [info@abielbiotech.com](mailto:info@abielbiotech.com).