

New Antibody for Gene Expression Research!!!

*Antibody against monoubiquitinated histone
H2B*

Antibody

Mouse Monoclonal Antibody, IgG2a isotype

Precise binding epitope

**This antibody was raised against a branched peptide
corresponding
to the conjugation site of ubiquitin on human histone H2B**

Tested Application

Western Blot, Chromatin Immunoprecipitation

Species Reactivity

Human, Mouse, Drosophila and Arabidopsis

Contact us:

5800 Royalmount, Montréal, Québec, Canada H4P 1K5
Tel.: (514) 940-MABS • Toll Free: (888) 940-MABS • Fax:
(514) 733-1212

www.medimabs.com info@medimabs.com