


Immunohistochemistry Reagents



Quality Antibodies for Quality Research

SouthernBiotech

www.SouthernBiotech.com info@SouthernBiotech.com



Immunohistochemistry Reagents

Fluoromount-G®

Mounting Medium	Cat. No.	Qty
Fluoromount-G®	0100-01	25 mL
DAPI Fluoromount-G®	0100-20	20 mL

Fluoromount-G® is a water-soluble, non-fluorescing compound for slides mounted after a staining procedure having an aqueous final step. It provides a semi-permanent seal for long-term storage of slide preparations.

Based on data from published papers, Fluoromount-G® mounting media has been shown to be compatible with the following fluorochromes and reagents

AEC	BODIPY® 493/503	CY3.5	eFluor® 615	Fluoro-Jade® C	MitoTracker® Orange	Propidium Iodide (PI)
Alexa Fluor® 350	BODIPY® FL	CY5	EGFP	Fluoro-Ruby®	MitoTracker® Red	PROTEOSTAT®
Alexa Fluor® 405	BODIPY® TR Glibenclamide	CY5.5	Eosin	FM® 4-64	MitoTracker® Red CMXRos	Qtracker® 655
Alexa Fluor® 488	Brilliant Violet 421™	DAB	ER-Tracker™ Green	GFP	mOrange	RFP
Alexa Fluor® 532	CdSe-ZnS Quantum Dots	DAPI	ER-Tracker™ Red	HCS LipidTOX™ Deep Red	mRFP	Rhodamine (TRITC)
Alexa Fluor® 546	CellMask™ Orange	DiD	Ethidium Bromide	HCS LipidTOX™ Green	MTT	Rhodamine 6G
Alexa Fluor® 555	CellROX® Green	DiI	Ethidium Monoazide (EMA)	Hematoxylin	NBD C6-Ceramide	Rhodamine B
Alexa Fluor® 568	CellTracker™ Blue	DiO	Evans Blue	Hoechst 33258	NBT	Rhodamine Red™-X
Alexa Fluor® 594	CellTracker™ CM-Dil	DRAQ5®	EYFP	Hoechst 33342	NeuroTrace® 435/455 Blue	Rhodamine Red™-X
Alexa Fluor® 633	CellTracker™ Green	DsRed	Fast Blue	Hoechst 34580	NeuroTrace® 640/660 Deep-Red	Sudan Black B
Alexa Fluor® 635	CellTracker™ Orange	dTomato	Fast Red	H-Ruby	Nile Red	SYTO™ Green
Alexa Fluor® 647	CF®488A	Duolink® Red	FD Rapid GolgiStain™	Image-iT® FX Signal Enhancer	NorthernLights™ 493	TACS® 2 TdT Blue
Alexa Fluor® 660	CF®568	DyLight® 350	FLIVO®	Lissamine™ Rhodamine B	NorthernLights™ 557	Texas Red®
Alexa Fluor® 680	CF®660C	DyLight® 405	FluoProbes® 488	LysoTracker® Green	Octadecyl Rhodamine B	Thioflavin S
Alizarin Red S	Chlorin e6 (Ce6)	DyLight® 488	FluoProbes® 547H	LysoTracker® Red	Oil Red O	TMR Red
AMCA	Chromoem™ 642	DyLight® 550	Fluorescein (FITC)	Marina Blue®	OPF	TO-PRO®-3
AmCyan	Coumarin	DyLight® 594	Fluoro-Emerald	mCherry	Oregon Green® 488	TSA™ CY3
Amylo-Glo®	CY2	DyLight® 633	Fluoro-Gold	MitoSOX® Red	Pacific Blue™	X-gal
APC	CY2	DyLight® 650	Fluoro-Jade®	MitoTracker® Deep Red	Phycocyanin (PE)	YFP
BCIP/NBT	CY3	eFluor® 450	Fluoro-Jade® B	MitoTracker® Green	PE/CY	

Microscope Slides

	Cat. No.	Qty
Gelatin Subbed Slides	SLD01-BX	72 Slides/Box
	SLD01-CS	10 Boxes/Case

SouthernBiotech's Gelatin Subbed Slides have been carefully cleaned and subbed with a solution of 0.1% gelatin and chromium potassium sulfate. The coated slides measure 25 mm x 75 mm with a frosted end for marking.

Corporate Offices

160 Oxmoor Blvd
Birmingham, AL 35209
USA

Mailing Address

PO Box 26221
Birmingham, AL 35260
USA

Tel 205.945.1774

Fax 205.945.8768

Email info@SouthernBiotech.com


Technical Support technical@SouthernBiotech.com

Website www.SouthernBiotech.com

SouthernBiotech is ISO 9001:2015 Certified

SouthernBiotech

www.SouthernBiotech.com info@SouthernBiotech.com



Fluoromount-G® is a registered trademark of Southern Biotechnology Associates, Inc.
Alexa Fluor®, BODIPY®, CellROX®, eFluor®, DyLight®, FM®, Image-iT®, LysoTracker®, Marina Blue®, MitoSOX®, MitoTracker®, Oregon Green®, NeuroTrace®, Qtracker®, SYTOX®, Texas Red®, and are TO-PRO® registered trademarks of Thermo Fisher Scientific, Inc. and its subsidiaries
Amylo-Glo®, Fluoro-Jade®, and Fluoro-Ruby® are registered trademarks of Histo-Chem, Inc.
CellMask™, CellTracker™, ER-Tracker™, LipidTOX™, Pacific Blue™, and Rhodamine Red™ are trademarks of Thermo Fisher Scientific, Inc. and its subsidiaries
FLIVO® is a registered trademark of Seed Research and Development, LLC
Duolink® is a registered trademark of Sigma-Aldrich International GmbH
FD Rapid GolgiStain™ is a trademark of FD NeuroTechnologies, Inc.
PROTEOSTAT® is a registered trademark of Enzo Life Sciences, Inc.
Lissamine™ is a trademark of Imperial Chemical Industries, Inc.
FluoProbes® is a registered trademark of INTERCHIM S.A.
Brilliant Violet 421™ is a trademark of Sirigen Group Ltd.
DRAQ5® is a registered trademark of Sigma-Aldrich International GmbH
NorthernLights™ is a trademark of B&D Systems, Inc.
TACS® is a registered trademark of Trevigen, Inc.
Cy® is a registered trademark of GE Healthcare
CF® is a registered trademark of Biotium, Inc.
Chromoem™ is a trademark of Active Motif
TSA™ is a trademark of PerkinElmer, Inc.
MM-IHC-02