

Horse-side results for Horse-side decisions



Intuitive blood sampling

from syringe or blood tube




Numerical results

for objective assessment & monitoring



No lost results

with last result saved by reader



Early identification of systemic inflammation and infection is challenging, as the clinical signs can be subtle. SAA is an indispensable tool to aid equine veterinarians in this challenge, since levels increase with acute inflammation even more rapidly and dramatically than elevations in fibrinogen or WBC count. SAA starts dropping as soon as disease begins to resolve, allowing objective monitoring of clinical condition and treatment efficacy.

Multiple tests, one handheld reader - now with RECHARGEABLE battery!

- ✓ Ideal size with convenient storage case
- ✓ Calibration for accuracy & consistency
- ✓ Will not lose your results
- ✓ Automatically reads test results at 10 minutes

