



InVivoMab vs. *InVivoPlus*

	<u><i>InVivoMab</i></u>	<u><i>InVivoPlus</i></u>
purity level	> 95%	> 95%
protein integrity	✓ (verified via SDS-PAGE)	✓ (verified via SDS-PAGE)
endotoxin concentration	< 2EU/mg	< 1EU/mg
validated by WB, FC, or ELISA		✓
protein aggregates validated at ≤ 5%		✓
tested for murine pathogens		✓
azide and carrier protein free	✓	✓
available in bulk quantities	✓	✓
suited for <i>in vivo</i> studies	✓	✓
catalog number	starts with BE	starts with BP

With more than 3000 product citations in quality peer-reviewed journals Bio X Cell is a trusted supplier of antibodies for academic, private, and non-profit institutions. These publications provide the best possible testimony for the exceptional quality of our products. Up-to-date reference lists for each of our products are available on our website. For more information on how we can help you achieve your research goals, please browse our website and contact us with any questions.



Bulk Antibodies for *in vivo* Research

1.866.787.344 | bxcell.com | customerservice@bxcell.com