



MitoFlamma[®] (live)



Mitochondria-selective fluorescent dye for living cell analysis

- Staining mitochondria regardless of mitochondrial potential
- Sustainable during the fixation step
- Sustainable after permeabilization steps.





Fluorescence intensity of MitoFlamma[®] Green (left) and Deep Red (right)

MitoFlamma [®]	λ_{Ex} (nm)	λ_{Em} (nm)
● MitoFlamma [®] Green	498	583
● MitoFlamma [®] Deep Red	630	660



MitoFlamma[®] Green





MitoFlamma[®] Deep Red



MitoFlamma[®] (live)

Mitochondria-selective fluorescent dye for living cell analysis



Cytotoxicity comparison of MitoFlamma[®] and other mitochondrial probe

Custom Service: BioActs provides custom antibody/protein fluorescence labeling services with a wide selection of reactive dyes and proper labeling technology

For ordering custom service, complete custom labeling service request form and submit via email order@bioacts.com or by fax to +82-32-818-8206

For questions, please contact us at support@bioacts.com

Headquarter
ADDRESS: 9, Cheongneung-daero 595beon-gil
Namdong-gu, Incheon 21666, Korea
WEBSITE: www.bioacts.com
PHONE: +82-32-818-9100
FAX: +82-32-818-8206

MAILS
ORDER SUPPORT: order@bioacts.com
CUSTOMER SUPPORT: support@bioacts.com
B2B/BULK ORDER SUPPORT: ivd@bioacts.com

www.bioacts.com