



SMALL MOLECULES

Antigen/Antibody

- Product List
- Vitamins, Hormones, Antibiotics, DOA Therapeutic Drugs, Mycotoxins, etc.



Small Molecules

Vitamins

Vitamin A
 Vitamin B1, B2, B3, B7
 Vitamin B7/Biotin
 Vitamin B9/Folic Acid
 Vitamin B12
 Vitamin C
 25-OH Vitamin D
 1,25 (OH)₂ Vitamin D3
 Vitamin E
 Vitamin K1

Hormones

19-Norandrosterone
 Bisphenol a
 Estradiol
 Hydrocortisone
 Methyltestosterone
 Nonylphenol
 Tetrabromobisphenol A
 Polybrominated diphenyl ether
 Ractopamine
 Testosterone

Antibiotics

Sulfadiazine Avermectin
 Sulfaquinoxaline Tylosin
 Sulfamethazine Tilmicosin
 Sulfadimethoxine Vancomycin
 Trimethoprim Salinomycin
 Chloramphenicol Monensin
 Kanamycin Tetracycline
 Streptomycin Olaquinox
 Neomycin Sarafloxacin
 Erythrocin Ofloxacin
 Spectinomycin Enrofloxacin
 Fluoroquinolones

Drugs of Abuse

Acetamidrid
 Acetochlor
 Carbaryl
 Carbendazim
 Chlorothalonil
 Chlorpyrifos
 Clonazepam
 Diazepam
 Dihydrochlorothiazide
 Fenpropathrin
 Imidacloprid
 Nitrazepam
 Parathion methyl

Therapeutic Drugs

Carbamazepine
 Digoxin
 L-Norephedrine
 Phenytoin
 Theophylline
 Valproic Acid

Mycotoxins

Aflatoxin B1
 Aflatoxin M1
 Aflatoxin
 Ochratoxin
 Zearalenone
 Fumonisin
 T-2 toxin


Amino Acids

Adrenaline
 Agmatine
 Anthranilic Acid
 Amino-levulinic Acid
 Carnosine
 Diamino-Pentane
 Dopamine
 GABA
 Glycine
 Histamine
 Homocysteic Acid
 Homocysteine
 Kainic Acid
 Kynurenic Acid
 NO₂-Tyrosine

Small Molecule Antigen and Antibody

Haptens, also known as incomplete antigens, are small molecules that can bind to the antibody, but will usually not initiate an immune response alone. It only has immunoreactivity, no immunogenicity. However, it can elicit an immune response and produce the antibody when attached to a large carrier such as a protein. Most polysaccharides and all lipids are haptens. Some small substances (such as vitamins, antibiotics, drugs, etc.) with specific structure called simple hapten. When simple hapten gets into the body of allergies, it can combine with the tissue protein and become complete antigen which can cause hypersensitivity. In general, B lymphocytes recognize hapten epitopes and T lymphocytes recognize carrier epitopes.

Small Molecules Antigen/Antibody



Vitamin

Specificity	Catalog	Type	Host	Application	Format
Vitamin A	DPATB-H83255	PAb	Rabbit	ELISA	N/A
Vitamin B ₁	DPATB-H83275	PAb	Rat	ELISA, IHC-Fr	N/A
Vitamin B ₂	DPAB1793	PAb	Rat	N/A	N/A
Vitamin B ₃	DPATB-H83259	PAb	Rat	ELISA	N/A
Vitamin B ₅	DPABHZ12	PAb	Rabbit	ELISA	N/A
Vitamin B ₆	DPATB-H83257	PAb	Rat	ELISA, IHC-Fr	N/A
Vitamin B ₇ /Biotin	DMABT-54801MB	MAB	Mouse	ELISA, IHC, WB	N/A
Vitamin B ₇ /Biotin	DPAB21486	PAb	Goat	ELISA, WB	HRP
Vitamin B ₉ /Folic Acid	DMABT-54808MF	MAB	Mouse	IHC, ELISA, FC, IF, IP, WB	N/A
Vitamin B ₉ /Folic Acid	DPAB1713	PAb	Rat	IHC, ELISA, FC, IF, IP, WB	N/A
Vitamin B ₁₂	DMAB6300	MAB	Mouse	ELISA	N/A
Vitamin B ₁₂	DPATB-H83238	PAb	Rabbit	ELISA	N/A
Vitamin C	DPATB-H83260	PAb	Rat	EM, ELISA	N/A
25-OH Vitamin D	DMAB2914	MAB	Mouse	ELISA	N/A
1,25 (OH) ₂ Vitamin D ₃	DMAB2913	MAB	Mouse	ELISA	N/A
Vitamin E	DPAB-DC3974	PAb	Rabbit	ELISA, IHC	N/A
Vitamin K ₁	DPABHZ13	PAb	Rabbit	ELISA, WB, IHC, IF/ICC	N/A

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Vitamin B ₇ /Biotin	HMABPY074	MAB	Mouse	ELISA, LFIA	Purified
Vitamin B ₉ /Folic Acid	HMABPY072	MAB	Mouse	ELISA, LFIA	Purified
Vitamin B ₁₂	HMABPY073	MAB	Mouse	ELISA, LFIA	Purified

Antigen

Catalog	Product Name	Application	Nature
DAG3411	Thiamine [G-BSA]	ELISA, LFIA	Synthetic
DAG3399	Riboflavin [BSA]	ELISA, LFIA	Synthetic
DAG3371	Nicotinamide [BSA]	ELISA, LFIA	Synthetic
DAG3394	Pyridoxal [BSA]	ELISA, LFIA	Synthetic
DAG3197	Biotin [BSA]	ELISA, LFIA	Synthetic
DAG2991	Folic Acid [BSA]	ELISA, LFIA	Synthetic
DAG3037	Vitamin B ₁₂ [BSA]	ELISA, LFIA	Synthetic
DAG03219	25-OH Vitamin D ₃ -BSA Conjugated	ELISA, LFIA	Synthetic
DAG3269	Alpha-Tocopherol [BSA]	ELISA, LFIA	Synthetic

Autoimmunity / Anti-Nuclear

Antibody

Specificity	Catalog	Type	Host	Application	Format
Centromere B	CABT-B8412	MAB	Mouse	WB, IP, IF	N/A
Centromere B	DPABH-06577	PAb	Mouse	WB, ICC/IF	N/A
dsDNA	DPATB-H82289	PAb	Sheep	IP, ISH	N/A
dsDNA	DMABT-PY6023	MAB	Mouse	Dot, ICC, IFA, ELISA	N/A
Histone H1	DPABH-01943	PAb	Rabbit	WB, ELISA, IHC-P, ICC/IF	N/A
Histone H2A	DMAB7957	MAB	Mouse	ELISA, Dot	N/A
Histone H2A	CABT-LM00282	PAb	Rabbit	WB, IP, IHC-P	N/A
Histone H2B	DPABH-12995	PAb	Rabbit	WB, IHC-P, ChIP, ICC/IF, IP	N/A
Histone H3	DMAB7958	MAB	Mouse	ELISA, Dot	N/A
Histone H3	DPAB-DC4361	PAb	Rabbit	ChIP, IF, ELISA, Dot	N/A
Histone H4	DCABH-8508	MAB	Mouse	FC, IHC-P, IP, WB, ChIP, ICC/IF	N/A
Histone H4	DPABH-06057	PAb	Rabbit	ICC/IF, IHC-P, EM, ChIP, WB	N/A
Jo-1	DCABH-3815	MAB	Rabbit	WB, IHC-P, IP	N/A
Jo-1	DPABH-00625	PAb	Rabbit	WB, IHC-P	N/A
Sci-70	DCABH-2459	MAB	Rabbit	WB, ICC/IF	N/A
Sci-70	DPAB-DC3092	PAb	Mouse	WB, ELISA	N/A
Sm	DMABT-Z60299	MAB	Mouse	EM, ICC/IF, WB, IHC-P	N/A
SSA	DCABH-4115	MAB	Mouse	WB, ELISA, IHC-P	N/A
SSA	DPAB-DC2966	PAb	Mouse	WB, ELISA	N/A
SSB	DPAB-DC2967	PAb	Mouse	WB, ELISA	N/A

Positive Control Antibody

Specificity	Catalog	Type	Host	Application	Format
dsDNA	DPAB0553	PAb	Human	ELISA	Purified
ssDNA	DPABHZ01	PAb	Human	ELISA	Purified
Centromere	DPABHZ10	PAb	Human	ELISA	Purified
Histones	DPABHZ09	PAb	Human	ELISA	Purified
Jo-1	DPABHZ07	PAb	Human	ELISA	Purified
Ribosomal P	DPABHZ08	PAb	Human	ELISA	Purified
RNP	DPABHZ02	PAb	Human	ELISA	Purified
Scl-70	DPABHZ06	PAb	Human	ELISA	Purified
Sm	DPABHZ03	PAb	Human	ELISA	Purified
SS-A/Ro	DPABHZ04	PAb	Human	ELISA	Purified
SS-B/La	DPABHZ05	PAb	Human	ELISA	Purified

Hormones

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Bisphenol a	HMABPY029	MAb	Mouse	ELISA, LFIA	Purified
Tetrabromobisphenol A	HMABPY030	MAb	Mouse	ELISA, LFIA	Purified
Polybrominated diphenyl ether	HMABPY031	MAb	Mouse	ELISA, LFIA	Purified
Nonylphenol	HMABPY032	MAb	Mouse	ELISA, LFIA	Purified
Ractopamine	HMABPY001	MAb	Mouse	ELISA, LFIA	Purified
Salbutamol	HMABPY002	MAb	Mouse	ELISA, LFIA	Purified
Clenbuterol	HMABPY003	MAb	Mouse	ELISA, LFIA	Purified
Clenbuterol residues	HMABPY004	MAb	Mouse	ELISA, LFIA	Purified
19-Norandrosterone	HMABPY019	MAb	Mouse	ELISA, LFIA	Purified
Methyltestosterone	HMABPY020	MAb	Mouse	ELISA, LFIA	Purified
Testosterone	HMABPY021	MAb	Mouse	ELISA, LFIA	Purified
Medroxyprogesterone	HMABPY022	MAb	Mouse	ELISA, LFIA	Purified
Progesterone	HMABPY023	MAb	Mouse	ELISA, LFIA	Purified
Hexadecadrol	HMABPY024	MAb	Mouse	ELISA, LFIA	Purified
Diethylstilbestrol	HMABPY025	MAb	Mouse	ELISA, LFIA	Purified
Hexoestrolum	HMABPY026	MAb	Mouse	ELISA, LFIA	Purified
Hydrocortisone	HMABPY027	MAb	Mouse	ELISA, LFIA	Purified
Estradiol	HMABPY028	MAb	Mouse	ELISA, LFIA	Purified

Antigen

Catalog	Product Name	DAG1089	Dienestrol(p-OH) [HRP]
DAG3411	Thiamine [BSA]	DAG1086	Dexamethasone(21) [HRP]
DAG1000	17a-Testosterone(3) [HRP]	DAG1090	Diethylstilbestrol(p-OH) [HRP]
DAG1001	17a-Trenbolone(3) [HRP]	DAG1092	Dihydrotestosterone [HRP]
DAG1022	Androstenedione(3) [HRP]	DAG1179	Ethinylestradiol(6) [HRP]
DAG1002	17b-Boldenone(3) [HRP]	DAG1181	17a-Estradiol(3) [HRP]
DAG1003	17b-19-Nortestosterone(3) [HRP]	DAG1185	Estriol(6) [HRP]
DAG1054	Chlormadinone acetate(3) [HRP]	DAG1195	Flumethasone [HRP]
DAG1056	Cimaterol [HRP]	DAG034	Clenbuterol (CLEN) [BSA]
DAG1072	Cortisol(4) [HRP]	DAG033	Ractopamine (RAC) [BSA]
DAG1073	Cortisone [HRP]	DAG3263	Abcsic Acid [BSA]
DAG1081	DHEA(17) [HRP]		

Antibiotics

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Sulfadiazine	HMABPY035	MAb	Mouse	ELISA, LFIA	Purified
3-Methyl-quinosaline-2 -carboxylicacid	HMABPY062	MAb	Mouse	ELISA, LFIA	Purified
Albendazole	HMABPY057	MAb	Mouse	ELISA, LFIA	Purified
Ampicillin	HMABPY068	MAb	Mouse	ELISA, LFIA	Purified
Avermectin	HMABPY051	MAb	Mouse	ELISA, LFIA	Purified
Cefalexin	HMABPY071	MAb	Mouse	ELISA, LFIA	Purified
Cefquinome	HMABPY070	MAb	Mouse	ELISA, LFIA	Purified
Ceftiofur	HMABPY069	MAb	Mouse	ELISA, LFIA	Purified
Chloramphenicol	HMABPY040	MAb	Mouse	ELISA, LFIA	Purified
Ciprofloxacin	HMABPY064	MAb	Mouse	ELISA, LFIA	Purified
Enrofloxacin	HMABPY067	MAb	Mouse	ELISA, LFIA	Purified
Erythrocin	HMABPY049	MAb	Mouse	ELISA, LFIA	Purified
Florfenicol	HMABPY042	MAb	Mouse	ELISA, LFIA	Purified
Fluoroquinolones	HMABPY082	MAb	Mouse	ELISA, LFIA	Purified
Gentamicin	HMABPY043	MAb	Mouse	ELISA, LFIA	Purified
Kanamycin	HMABPY044	MAb	Mouse	ELISA, LFIA	Purified
Leucomalachite green	HMABPY094	MAb	Mouse	ELISA, LFIA	Purified
Lincomycin	HMABPY047	MAb	Mouse	ELISA, LFIA	Purified

Malachite green	HMABPY093	MAb	Mouse	ELISA, LFIA	Purified
Metronidazole	HMABPY059	MAb	Mouse	ELISA, LFIA	Purified
Monensin	HMABPY056	MAb	Mouse	ELISA, LFIA	Purified
Natamycin	HMABPY058	MAb	Mouse	ELISA, LFIA	Purified
Neomycin	HMABPY046	MAb	Mouse	ELISA, LFIA	Purified
Ofloxacin	HMABPY066	MAb	Mouse	ELISA, LFIA	Purified
Olaquinox	HMABPY061	MAb	Mouse	ELISA, LFIA	Purified
Pentachlorophenol	HMABPY095	MAb	Mouse	ELISA, LFIA	Purified
Pirlimycin	HMABPY048	MAb	Mouse	ELISA, LFIA	Purified
Quinolones	HMABPY063	MAb	Mouse	ELISA, LFIA	Purified
Salinomycin	HMABPY055	MAb	Mouse	ELISA, LFIA	Purified
Sarafloxacin	HMABPY065	MAb	Mouse	ELISA, LFIA	Purified
Spectinomycin	HMABPY050	MAb	Mouse	ELISA, LFIA	Purified
Streptomycin	HMABPY045	MAb	Mouse	ELISA, LFIA	Purified
Sulfadimethoxine	HMABPY038	MAb	Mouse	ELISA, LFIA	Purified
Sulfamethazine	HMABPY037	MAb	Mouse	ELISA, LFIA	Purified
Sulfanilamides	HMABPY083	MAb	Mouse	ELISA, LFIA	Purified
Sulfaquinoxaline	HMABPY036	MAb	Mouse	ELISA, LFIA	Purified
Tetracycline	HMABPY060	MAb	Mouse	ELISA, LFIA	Purified
Thiamphenicol drug Residue	HMABPY041	MAb	Mouse	ELISA, LFIA	Purified
Tilmicosin	HMABPY053	MAb	Mouse	ELISA, LFIA	Purified
Trimethoprim	HMABPY039	MAb	Mouse	ELISA, LFIA	Purified
Tylosin	HMABPY052	MAb	Mouse	ELISA, LFIA	Purified
Vancomycin	HMABPY054	MAb	Mouse	ELISA, LFIA	Purified

Antigen

Catalog	Product Name	Catalog	Product Name
DAG1018	Amikacin [HRP]	DAG1193	Florfenicol Amine [HRP]
DAG4461	Amoxicillin [BSA]	DAG1194	Flumequine [HRP]
DAG4477	Amoxicillin [KLH]	DAG1198	Anti-Gentamicin [HRP]
DAGB332	AMOX [HRP]	DAG2993	Gentamicin [BSA]
DAG1021	Ampicillin [HRP]	DAG4484	Gentamicin [KLH]
DAG4462	Ampicillin [BSA]	DAG4469	Hygromycin B [BSA]
DAG4478	Ampicillin [KLH]	DAG4485	Hygromycin B [KLH]
DAGB333	AOZ [HRP]	DAG1208	Kanamycin [HRP]
DAG1041	Apramycin [HRP]	DAG1213	Leucomalachite green [HRP]
DAG4463	Apramycin [BSA]	DAGB305	Malachite Green [HRP]
DAG4479	Apramycin [KLH]	DAG1246	Nalidixic Acid [HRP]
DAGB306	Avilamycin [HRP]	DAG1248	Neomycin [HRP]
DAG1043	Bacitracin [HRP]	DAG4470	Neomycin [BSA]
DAG4464	Bacitracin [BSA]	DAG4486	Neomycin [KLH]
DAG4480	Bacitracin [KLH]	DAG1011	4-Nitrophenyl AHD [HRP]
DAG1049	Carbadox [HRP]	DAGB309	4-NITROPHENYL-AOZ [HRP]
DAGB300	Olaquinox [HRP]	DAG1259	Norfloxacin [HRP]
DAG1052	Ceftiofur [HRP]	DAGB303	Ormetoprim [HRP]
DAGB302	Cephalexin [HRP]	DAG1267	Oxolinic Acid [HRP]
DAG1053	Chloramphenicol [HRP]	DAG1277	Pirlimycin [HRP]
DAG2974	Chloramphenicol [BSA]	DAG1285	Ronidazole [HRP]
DAG4436	Chloramphenicol [OVA]	DAG1291	Spiramycin [HRP]
DAG4467	Cloxacillin [BSA]	DAG1293	Streptomycin [HRP]
DAG4483	Cloxacillin [KLH]	DAG4472	Streptomycin [BSA]
DAGB335	Dicloxacillin [HRP]	DAG4488	Streptomycin [KLH]
DAGB308	Dinitrosalicylic Acid Hydrazide [HRP]	DAG1294	Sulfabenzamide v [HRP]
DAGB301	Erythromycin [HRP]	DAGB336	Sulfachloropyridazine [HRP]
DAGB307	Flavomycin [HRP]	DAG1297	Sulfadiazine [HRP]

Drugs of Abuse (DOA)

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Acetamidiprid	HMABPY090	MAb	Mouse	ELISA; LFIA	Purified
Acetochlor	HMABPY089	MAb	Mouse	ELISA; LFIA	Purified
Carbaryl	HMABPY088	MAb	Mouse	ELISA; LFIA	Purified
Carbendazim	HMABPY084	MAb	Mouse	ELISA; LFIA	Purified
Chlorothalonil	HMABPY085	MAb	Mouse	ELISA; LFIA	Purified
Chlorpyrifos	HMABPY087	MAb	Mouse	ELISA; LFIA	Purified
Clonazepam	HMABPY102	MAb	Mouse	ELISA; LFIA	Purified
Diazepam	HMABPY103	MAb	Mouse	ELISA; LFIA	Purified
Dihydrochlorothiazide	HMABPY104	MAb	Mouse	ELISA; LFIA	Purified
Fenprothrin	HMABPY092	MAb	Mouse	ELISA; LFIA	Purified
Imidacloprid	HMABPY091	MAb	Mouse	ELISA; LFIA	Purified
Nitrazepam	HMABPY101	MAb	Mouse	ELISA; LFIA	Purified
Parathion methyl	HMABPY086	MAb	Mouse	ELISA; LFIA	Purified

Antigen

Catalog	Product Name		
DAG1034	2-Amino-Benzimidazole Sulphone [HRP]	DAG1094	Dipyron [HRP]
DAG1015	Acepromazine [HRP]	DAG2983	EDDP [BTG]
DAG1161	Acepromazine [HRP]	DAG1159	Escitalopram [HRP]
DAG2957	Acetaminophen [BSA]	DAG2989	Ethyl Glucuronide [BSA]
DAG3275	Acetyl Salicylic Acid [BSA]	DAG1189	Ethyl glucuronide [HRP]
DAG3266	Acrolein [BSA]	DAG1190	Etodolac [HRP]
DAG1016	Albendazole [HRP]	DAG398	Fentanyl [BSA]
DAGB342	Aminoflunitrazepam [HRP]	DAG1122	Fentanyl [HRP]
DAGB330	Aminopyrine [HRP]	DAG1196	Flunixin(COOH) [HRP]
DAG409	Amphetamine (AMP) [BSA]	DAG1199	Halofuginone [HRP]
DAG1033	Benzimidazole [HRP]	DAG1164	Haloperidol [HRP]
DAG408	Benzodiazepine (BZO) [BSA]	DAG1167	Heroin [HRP]
DAG396	Benzoylcegonine [BSA]	DAG1203	Hydrocodone(6) [HRP]
DAG413	Benzoylcegonine [HRP]	DAG3353	Ibuprofen [BSA]
DAG394	Buprenorphine [BSA]	DAG1327	Imipramine [HRP]
DAG393	Buprenorphine [BTG]	DAG1209	Kebuzone [HRP]
DAG419	Buprenorphine [HRP]	DAG400	Ketamine [BSA]
DAG3099	Cannabinoid THC [BSA]	DAG1211	Ketoprofen [HRP]
DAG1162	Carboxymethyl LSD [HRP]	DAG1125	LAAM [HRP]
DAG3107	Carisoprodol [BSA]	DAG1212	Lasalocid A [HRP]
DAG1157	Chlorpromazine [HRP]	DAG1216	Lorazepam(3) [HRP]
DAG1123	Codeine(6) [HRP]	DAG395	MDMA [BTG]
DAG417	Cotinine-3 [HRP]	DAG1223	Meloxicam [HRP]
DAG418	Cotinine-4 [HRP]	DAG1224	Meperidine [HRP]
DAG3100	C-THC gen [BSA]	DAG1225	Meprobamate [HRP]
DAG1078	Cyclopentobarbital (5) [HRP]	DAG1226	Mescaline [HRP]
DAG1312	D8-THC-COOH [HRP]	DAG2998	Methadone [BSA]
DAG1079	Dapsone [HRP]	DAG1124	Methadone(p) [HRP]
DAG1037	Decoquinat [HRP]	DAG422	Methamphetamine-N [BSA]
DAG2980	Delta 9 THC [BSA]	DAG1242	Monensin [HRP]
DAGB329	Desmethylnaproxen [HRP]	DAG420	Morphine [BSA]
DAG1168	Dextromethorphan [HRP]	DAG3009	MORPHINE [BSA]
DAG1026	Diclazuril [HRP]	DAG3140	Morphine [BTG]
DAG3299	Diclofenac [BSA]	DAG1243	Morphine [HRP]
DAG1091	Digoxigenin [HRP]	DAG1247	Naltrexone [HRP]
		DAG3370	Naproxen [BSA]

DAG1249	Nicarbazin [HRP]
DAG1251	Nitrazepam(7) [HRP]
DAG3373	NO-Acetyl Salicylic Acid [BSA]
DAG1256	Norbuprenorphine(3) [HRP]
DAG1257	Nordiazepam [HRP]
DAG3106	Oxazepam [BSA]
DAG1265	Oxazepam(3) [HRP]
DAG3014	Oxycodone [BSA]
DAG1268	Oxycodone(N) [HRP]
DAG1275	Oxyphenbutazone [HRP]
DAG3015	PCP [BSA]
DAG3016	PCP [BTG]
DAG1270	PCP(m) [HRP]
DAG3103	Phencyclidine
DAG3017	Phencyclidine (PCP) [BSA]
DAG1273	Phenobarbital(m/p) [HRP]
DAGB331	Phenylbutazone [HRP]
DAG397	Propoxyphene [BSA]
DAG3022	Propoxyphene [HRP]

DAG3400	Rotenone [BSA]
DAG1289	Salvinorin [HRP]
DAG3102	Secobarbital [BSA]
DAG1290	Secobarbital [HRP]
DAG3026	THC [BSA]
DAG423	THC Delta 8 [BSA]
DAG414	THC Delta 8 [HRP]
DAG403	THC Delta 8+9 [BTG]
DAG404	THC Delta 9 [BSA]
DAG1027	Thiabendazole [HRP]
DAG1315	Thiabendazole [HRP]
DAG1318	Tolfenamic Acid [HRP]
DAG1319	Toltrazuril [HRP]
DAG399	Tramadol [BSA]
DAG1321	Tramadol [HRP]
DAG1323	Trazodone [HRP]
DAG1326	Triclabendazole [HRP]
DAG1336	Zolpidem [HRP]

Therapeutic Drugs

Antibody

Specificity	Catalog	Type	Host	Application	Format
Carbamazepine	DMAB3100	MAb	Mouse	ELISA	Purified
Carbamazepine	DPAB3468	PAb	Sheep	ELISA	Purified
Carbamazepine	DMABT-54912MC	MAb	Mouse	ELISA	Purified
Carbamazepine	DPBT-68253SC	PAb	Sheep	ELISA	Purified
Digoxigenin	DPBT-68262SD	PAb	Sheep	ELISA	Purified
Digoxin	DMABT-Z60491	MAb	Mouse	ELISA	FITC-Conjugated
Digoxin	DMAb6162	MAb	Mouse	ELISA	Purified
Digoxin	DPATB-H83017	PAb	Goat	ELISA	Purified
Digoxin	DPATB-H83095	PAb	Duck	ELISA	Purified
Digoxin	DPATB-H83405	PAb	Rabbit	ELISA	Purified
Digoxin	DPABH-00043	PAb	Sheep	ELISA	Purified
Digoxin	DPBT-68263SD	PAb	Sheep	ELISA	Purified
Diln	DPAB3471	PAb	Sheep	ELISA	Purified
Human Theophylline	DMAB2217	MAb	Mouse	ELISA	Purified
L-Norephedrine	DPAB-DC4502	PAb	Rabbit	ELISA	Purified
Methotrexate	DPATB-H81809	PAb	Goat	ELISA	Purified
Methotrexate	DPABY-070	PAb	Sheep	ELISA	Purified
Mycophenolic acid	DPAB3637	PAb	Sheep	ELISA	Purified
N-acetylprocainamide	DMABT-Z59839	MAb	Mouse	ELISA	Purified
Nitrendipine	DPAB-DC4570	PAb	Rabbit	ELISA	Purified
Nitrendipine	DPATB-H82326	PAb	Rabbit	ELISA	Purified
Norpropoxyphene	DPABY-931	PAb	Sheep	ELISA	Purified
Phenytoin	DMABT-Z60466	MAb	Mouse	ELISA	Purified
Phenytoin	DPATB-H82984	PAb	Rabbit	ELISA	Purified
Phenytoin	DPATB-H83000	PAb	Duck	ELISA	Purified
Phenytoin	DPAB4185	PAb	Sheep	ELISA	Purified
Phenytoin	DMABT-54918MP	MAb	Mouse	ELISA	Purified
Procainamide	DMABT-Z60407	MAb	Mouse	ELISA	Purified
Propranolol	DPABY-087	PAb	Sheep	ELISA	Purified
Quinidine	DPATB-H81934	PAb	Sheep	ELISA	Purified
Quinidine	DMABT-Z59931	MAb	Mouse	ELISA	Purified

Theophylline	DPATB-H82190	PAb	Duck	ELISA	Purified
Theophylline	DPABY-995	PAb	Sheep	ELISA	Purified
Theophylline	DMAB4380	MAB	Mouse	ELISA	Purified
Theophylline	DPAB27247	PAb	Rabbit	ELISA	Purified
Theophylline	DPAB3486	PAb	Sheep	ELISA	Purified
Theophylline	DPAB3487	PAb	Rabbit	ELISA	Purified
Valproic acid	DPAB-DC4866	PAb	Rabbit	ELISA	Purified
Valproic Acid	DMAB5330	MAB	Mouse	ELISA	Purified
Warfarin	DPABY-109	PAb	Rabbit	ELISA	Purified
Warfarin	DPAB28029	PAb	Sheep	ELISA	Purified

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Ribavirin	HMABPY033	MAB	Mouse	ELISA; LFIA	Purified
Amantadine	HMABPY034	MAB	Mouse	ELISA; LFIA	Purified

Antigen

Catalog	Product Name	Application	Nature
DAG1166	Carbamazepine [HRP]	ELISA, LFA	synthetic
DAG3043	Digoxin [BSA]	ELISA, LFA	synthetic
DAG172	Digoxin [FITC]	ELISA, LFA	synthetic
DAG3044	Digoxin [HRP]	ELISA, LFA	synthetic
DAG3365	L-Norephedrine [BSA]	ELISA, LFA	synthetic
DAG1276	Phenytoin [HRP]	ELISA, LFA	synthetic
DAG426	Theophylline [BSA]	ELISA, LFA	synthetic
DAG3030	Theophylline [HRP]	ELISA, LFA	synthetic
DAG3422	Valproic Acid [BSA]	ELISA, LFA	synthetic

Mycotoxins

Antibody

Specificity	Catalog	Type	Host	Application	Format
Aflatoxin B1/B2/Mycotoxins	DCABY-838	MAb	Mouse	ELISA	Purified
Aflatoxin B1/Mycotoxins	DMABT-Z59171	MAb	Mouse	ELISA	Purified
Aflatoxin B1-DNA Adduct/Mycotoxins	DMAB8913	MAb	Mouse	ELISA, IF, IHC	Purified
Aflatoxin/Mycotoxins	DMAB8910	MAb	Mouse	ELISA	Purified
Aflatoxin/Mycotoxins	DMABT-Z59170	MAb	Rat	ELISA	Purified
B-Zearalanol/Mycotoxins	DPAB-DC4032	PAb	Sheep	EIA	Purified
Deoxynivalenol/Mycotoxins	DPAB31452	PAb	Sheep	N/A	Purified
DON/Mycotoxins	DPAB4250	PAb	Rabbit	N/A	Purified
Fumonisin/Mycotoxins	DMAB8912	MAb	Mouse	ELISA	Purified
Fumonisin/Mycotoxins	DPAB31456	PAb	Rabbit	N/A	Purified
HT-2 toxin/Mycotoxins	DMAB21306	MAb	Mouse	N/A	Purified
Ochratoxin A/Mycotoxins	DPATB-H82292	PAb	Rabbit	ELISA	Purified
Ochratoxin A/Mycotoxins	DMAB8909	MAb	Mouse	ELISA	Purified
Patulin/Mycotoxins	DPAB-CL010	PAb	Rabbit	ELISA	Purified
ZEA/Mycotoxins	DMAB8911	MAb	Mouse	ELISA	Purified
Zearalanol/Mycotoxins	DPAB28030	PAb	Rabbit	IA, ELISA, WB, Dot, IHC	Purified

RHA™ mAb

Specificity	Catalog	Type	Host	Application	Format
Aflatoxin B1	HMABPY009	MAb	Mouse	ELISA, LFIA	Purified
Aflatoxin	HMABPY011	MAb	Mouse	ELISA, LFIA	Purified
Aflatoxin M1	HMABPY010	MAb	Mouse	ELISA, LFIA	Purified
Deoxynivalenol	HMABPY016	MAb	Mouse	ELISA, LFIA	Purified
Fumonisin	HMABPY014	MAb	Mouse	ELISA, LFIA	Purified
Ochratoxin	HMABPY012	MAb	Mouse	ELISA, LFIA	Purified
Sterigmatocystin	HMABPY017	MAb	Mouse	ELISA, LFIA	Purified
T-2 toxin	HMABPY015	MAb	Mouse	ELISA, LFIA	Purified
Zearalenone	HMABPY013	MAb	Mouse	ELISA, LFIA	Purified

Antigen

Catalog	Product Name	Application	Nature
DAGHZ05	Aflatoxin [BSA]	ELISA, LFA	synthetic
DAGHZ03	Aflatoxin B1 [BSA]	ELISA, LFA	synthetic
DAGHZ04	Aflatoxin M1 [BSA]	ELISA, LFA	synthetic
DAGHZ08	Deoxynivalenol [BSA]	ELISA, LFA	synthetic
DAGHZ06	Fumonisin [BSA]	ELISA, LFA	synthetic
DAG4437	Ochratoxin A [HRP]	ELISA, LFA	synthetic
DAGHZ09	Sterigmatocystin [BSA]	ELISA, LFA	synthetic
DAGHZ07	T-2 toxin [BSA]	ELISA, LFA	synthetic
DAG1333	Zearalanone(7) [HRP]	ELISA, LFA	synthetic

