




Isotype Control

APPLICATIONS:

- research use only
- Western blot
- Immunoassays

INFORMATION:

- InVivo BioTech Services GmbH

is certified according to -
ISO 9001:2015 - ISO
13485:2016

- Bulk production up to 1 kg

InVivo BioTech Services GmbH is a leading European biomanufacturer with ISO 9001 and ISO 13485 certification. We have been a Bruker company since 2017.

For more than 25 years, InVivo has been committed to supporting academics and industries with our competence, expertise, and reliable products.

We're here to support you and give you the confidence to bring your ideas to life with greater efficiency.

Isotype controls are essential reagents used in flow cytometry and immunoassays to differentiate non-specific background signals from true antigen-specific bindings. By matching the antibody class, subclass, and labeling, they ensure the accuracy and reliability of your experimental results.



Product Name	Cat. No.	Isotype	Aliquots
Mouse IgG2b k Isotype Control (MPC-11), InVivoPure (≤ 2 EU/mg)	<u>PSZ985P</u>	Mouse IgG2b k	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG2b k Isotype Control (MPC-11), InVivoPure+ (≤ 1 EU/mg)	<u>PSZ985P+</u>	Mouse IgG2b k	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG1 κ Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3421P</u>	Mouse IgG1 κ	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG1 κ Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3421P+</u>	Mouse IgG1 κ	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG2a κ Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3399P</u>	Mouse IgG2a κ	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG2a κ Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3399P+</u>	Mouse IgG2a κ	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG2b k Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3456P</u>	Mouse IgG2b k	1 mg 5 mg 25 mg 50 mg 100 mg
Mouse IgG2b k Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3456P+</u>	Mouse IgG2b k	1 mg 5 mg 25 mg 50 mg 100 mg
Rat IgG2b Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3616P</u>	Rat IgG2b	1 mg 5 mg 25 mg 50 mg 100 mg
Rat IgG2b Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3616P+</u>	Rat IgG2b	1 mg 5 mg 25 mg 50 mg 100 mg
Rat IgG2a Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3617P</u>	Rat IgG2a	1 mg 5 mg 25 mg 50 mg 100 mg
Rat IgG2a Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3617P+</u>	Rat IgG2a	1 mg 5 mg 25 mg 50 mg 100 mg



Product Name	Cat. No.	Isotype	Aliquots
Rat IgG1 Isotype Control, InVivoPure (≤ 2 EU/mg)	<u>AK3619P</u>	Rat IgG1	1 mg 5 mg 25 mg 50 mg 100 mg
Rat IgG1 Isotype Control, InVivoPure+ (≤ 1 EU/mg)	<u>AK3619P+</u>	Rat IgG1	1 mg 5 mg 25 mg 50 mg 100 mg



IVS CS-0136 Rev. 02